

AmpliAq®*Applied Biosystems***DNA Polymerase**

**with GeneAmp® 10X PCR Buffer
or GeneAmp® 10X PCR Buffer II
& MgCl₂ Solution**

**FOR RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC PROCEDURES**
See notice to purchaser

Enzyme	Enzyme with Buffer	
	GeneAmp 10X PCR Buffer	GeneAmp 10X PCR Buffer II & MgCl ₂ Solution
AmpliAq DNA Polymerase 250 Units (U), 5 U/μL, 200 rxns	N808-0160	N808-0161
AmpliAq DNA Polymerase 1000 Units (U), 5 U/μL, 800 rxns	N808-0171	N808-0172
Six Paq, AmpliAq DNA Polymerase 6 x 250 U, 5 U/μL, 1200 rxns	N808-0166	N808-0167
Twelve Paq, AmpliAq DNA Polymerase 12 x 250 U, 5 U/μL, 2400 rxns	N808-0152	N808-0153
AmpliAq DNA Polymerase, 5000 U 5 x 1000 U, 5 U/μL, 4000 rxns	N808-0155	N808-0156
AmpliAq DNA Polymerase, 25,000 U 25 x 1000 U, 5 U/μL, 20,000 rxns	N808-0185	N808-0186

TABLE OF CONTENTS

Notice to Purchaser 2
Storage and Stability..... 2
Introduction..... 3
List of Components..... 3
Enzyme Characteristics 4
Materials Required But Not Supplied..... 5
Protocols for DNA Amplification..... 5
1.0 General Advice..... 5
2.0 Reaction Mix Preparation..... 6
3.0 Amplification of DNA Templates..... 7
4.0 Temperature Cycling and Cycle Optimization..... 8
Performance Characteristics 9
Troubleshooting 11
References 13
Ordering Information..... 16

STORAGE AND STABILITY

Upon receipt, store the AmpliTaq DNA Polymerase and reagents at $-20^{\circ}\text{C} \pm 10^{\circ}\text{C}$ in a constant temperature freezer. If stored under the recommended conditions, the product will maintain performance through the control date printed on the label.

INTRODUCTION

Applied Biosystems AmpliTaq DNA Polymerase is a recombinant, thermostable, 94 kDa DNA polymerase encoded by a modified form of the *Thermus aquaticus* DNA polymerase gene which has been inserted into an *Escherichia coli* host.¹ This ultrapure, stable polymerase is designed for performing the Polymerase Chain Reaction (PCR) process and is available either by itself or as a component of the GeneAmp PCR Reagent Kit (Part No. N801-0055) and the GeneAmp PCR Core Reagents (Part No. N808-0009). For manual cycle sequencing of DNA, use the AmpliCycle™ Sequencing Kit (N808-0175) featuring AmpliTaq DNA Polymerase, CS.

The PCR process is a simple and powerful method,⁵ invented by K. Mullis, which allows amplification of DNA segments (including cDNA)⁶ *in vitro* through a succession of incubation steps at different temperatures.^{2,4,7,8,9} Typically, the double-stranded DNA is heat-denatured, the two primers complementary to the 3' boundaries of the target segment are annealed at low temperature and then extended at an intermediate temperature. One set of these three consecutive steps is referred to as a cycle. The PCR process is based on the repetition of this cycle and can amplify DNA segments by at least 10⁵ fold, and potentially as high as 10⁹ fold under the conditions described below.²

AmpliTaq DNA Polymerase gives optimal results under reagent conditions provided by GeneAmp 10X PCR Buffer and GeneAmp 10X PCR Buffer II. These buffers provide preferred pH and ionic strength for PCR amplification reactions (see below for exact formulations). The magnesium ion concentration required to achieve optimal PCR amplification is dependent on the specific set of primers and template used. GeneAmp 10X PCR Buffer contains 15 mM MgCl₂ which would give a 1.5 mM final concentration when diluted to 1X. The 25 mM MgCl₂ solution supplied with GeneAmp 10X PCR Buffer II can be used to adjust magnesium ion concentration for any set of primer-template pairs.

LIST OF COMPONENTS

Reagent	Volume	Description
AmpliTaq DNA Polymerase	50 µL	1 tube, 250 U of 5 U/µL AmpliTaq DNA Polymerase.
	— or — 200 µL	1 tube, 1000 U of 5 U/µL AmpliTaq DNA Polymerase.
GeneAmp 10X PCR Buffer (N808-0006) (N808-0129)	1.5 mL 6 x 1.5 mL	Tube(s) containing 100 mM Tris-HCl, pH 8.3, 500 mM KCl, 15 mM MgCl ₂ , 0.01% (w/v) gelatin.
GeneAmp 10X PCR Buffer II (N808-0010) (N808-0130) and 25 mM MgCl ₂ Solution	1.5 mL 6 x 1.5 mL and 1.5 mL 6 x 1.5 mL	Tube(s) containing 100 mM Tris-HCl, pH 8.3, 500 mM KCl. and Tube(s) containing 25 mM MgCl ₂ .
Product Insert	—	Containing AmpliTaq DNA Polymerase protocols.

Applied Biosystems

ENZYME CHARACTERISTICS

The key component is AmpliTaq DNA Polymerase (Deoxynucleosidetriphosphate: DNA Deoxynucleotidyltransferase, EC 2.7.7.7).

Concentration: 5 Units/ μ L

Unit Definition: One unit of enzyme is defined as the amount that will incorporate 10 nmoles of dNTPs into acid insoluble material per 30 minutes in a 10 minute incubation at 74°C under the analysis conditions below.

Analysis Conditions:¹ 25 mM TAPS (tris-[hydroxymethyl]-methyl-amino-propane-sulfonic acid, sodium salt), pH 9.3 (at room temperature); 50 mM KCl; 2 mM MgCl₂; 1 mM β -mercaptoethanol; 200 μ M each dATP, dGTP, dTTP; 100 μ M [α -³²P]-dCTP (0.05-0.1 Ci/mmol); salmon sperm DNA, activated by a modification of methods;¹⁰ mixed in a final volume of 50 μ L and incubated at 74°C for 10 minutes.

Storage Buffer: 20 mM Tris-HCl, pH 8.0 (at room temperature), 100 mM KCl, 0.1 mM EDTA (ethylenediaminetetraacetic acid), 1 mM DTT (dithiothreitol), 50% glycerol, 0.5% Tween® 20, 0.5% Nonidet® P 40.

Storage Temperature: Store AmpliTaq DNA Polymerase at -20°C, in a constant temperature freezer. If stored under proper conditions, the enzyme will remain active through the control date printed on the label.

Associated Activities: Endonuclease and exonuclease activities were not detectable after one hour incubation of 600 ng of supercoiled pBR322 (dam⁻, dcm⁻) or 600 ng of *Msp*I-digested pBR322 DNA, respectively, at 74°C, in the presence of 8 units of AmpliTaq DNA Polymerase. The enzyme has a fork-like-structure dependent, polymerization enhanced 5' to 3' nuclease activity but lacks a 3' to 5' exonuclease activity.^{3,11,12}

Applied Biosystems

MATERIALS REQUIRED BUT NOT SUPPLIED

Unless otherwise noted, the items listed are available from major laboratory suppliers (MLS) such as Baxter Scientific Products (McGaw Park, IL), Fisher Scientific (Pittsburgh, PA) or VWR Scientific (So. Plainfield, NJ).

Reagents

Autoclaved deionized, ultrafiltered
or glass distilled water, referred
to as "water" in this product insert.

Source

MLS

Equipment

Electrophoresis power supplyMLS

Gel electrophoresis equipmentMLS

GeneAmp® PCR Instrument System **Applied Biosystems**
(See page 16 of this
product insert).

Pipettors adjustableMLS

SuppliesPCR Reaction Tubes (see **Table 1** for tube/thermal cycler combination)GeneAmp Thin-Walled Reaction **Applied Biosystems,**
Tubes (0.5 mL, polypropylene) **(Part Nos. N801-0737,**
N801-0611, N801-0537).MicroAmp® Reaction Tubes**Applied Biosystems,**
(0.2 mL, polypropylene) **(Part Nos. N801-0533,**
N801-0540, N801-0612).Pipet tips with hydrophobic filtersMLS
used with adjustable pipettors**PROTOCOLS FOR DNA AMPLIFICATION****1.0 GENERAL ADVICE**

Due to the enormous amplification possible with the PCR process, small levels of DNA contamination, especially from previous PCR amplification reactions, samples with high DNA levels, and positive control templates, can result in product formation even in the absence of purposefully added template DNA.¹³ If possible, all reactions should be set up in an area separate from PCR product analysis. The use of dedicated or disposable vessels, solutions and pipettes (preferably positive displacement pipettes or tips with hydrophobic filters) for DNA preparation, reaction mixing, and sample analysis will minimize cross contamination.¹⁴ Use of the GeneAmp PCR Carry-over Prevention Kit (**Part No. N808-0068**) will prevent the reamplification of PCR products from previous amplifications.

Applied Biosystems

2.0 REACTION MIX PREPARATION

2.1 A Master Mix of reagents (water, buffer, dNTPs, primers and enzyme) for all samples can be prepared first, then aliquoted to individual tubes. Magnesium chloride and the template DNA are then added. Using such mixes will minimize reagent pipetting losses, increase accuracy, and reduce the number of reagent transfers. Perform amplifications in the Applied Biosystems PCR reaction tube appropriate for your GeneAmp PCR Instrument System (see table in **Section 4.1**, and the instrument manual). Applied Biosystems PCR reaction tubes provide the best heat transfer because of their uniform fit in the wells of the corresponding instrument.

2.2 Since DNA may stick to plastic and since nucleases are often found on surfaces, it may be preferable to use sterile, siliconized tubes and pipette tips.

2.3 Reaction Mix

	<u>Component</u>	<u>Addition</u>		<u>Final Concentration</u>
		<u>Order</u>	<u>Reaction</u>	
	Autoclaved deionized, ultrafiltered or glass distilled water	1	*	1X
Master Mix	10X PCR Buffer II (Part No. N808-0010)	2	5 μ L	10 mM Tris-HCl, pH 8.3, 50 mM KCl
	dATP	3	$\left\{ \begin{array}{l} 1 \mu\text{L} \\ 1 \mu\text{L} \\ 1 \mu\text{L} \\ 1 \mu\text{L} \end{array} \right.$	$\left. \begin{array}{l} 200 \mu\text{M} \\ 200 \mu\text{M} \\ 200 \mu\text{M} \\ 200 \mu\text{M} \end{array} \right\}$ or make daily as a mix of 1.25 mM of each dNTP; add 16 μ L per reaction
	dCTP 10 mM			
	dGTP each dNTP			
	dTTP			
	Primer 1	4	0.5-2.5 μ L	0.2-1.0 μ M
	Primer 2	5	0.5-2.5 μ L	0.2-1.0 μ M
AmpliAq DNA Polymerase	6	0.25 μ L	1.25 Units/50 μ L	
25 mM MgCl ₂ Solution	7	2-8 μ L	1.0-4.0 mM**	
Experimental template		*	500 ng/50 μ L	
Total Mix			50 μ L	

* Any combination of water and experimental template volumes can be used as long as the total volume of the reaction (including buffer, dNTP primers, enzyme and MgCl₂ solution) equals 50 μ L.

** The optimal magnesium chloride concentration may vary, depending on the primer and templates used and must be determined empirically. In most cases a final concentration of magnesium chloride in the range of 1.0 to 4.0 mM in the reaction mix will work well.

3.0 AMPLIFICATION OF DNA TEMPLATES

3.1 Gently mix (avoid generating bubbles) the AmpliTaq DNA Polymerase as well as other recently thawed reagents, then spin down in a microcentrifuge before pipetting. Pipette enzyme carefully and slowly; the viscosity of the 50% glycerol in the enzyme storage buffer can lead to errors. If possible use a positive displacement pipette. Using Master Mixes (see **Section 2.0**) will increase accuracy, reduce reagent loss on tips, and reduce tube-to-tube variability.

3.2 Except for amplification in the GeneAmp PCR System 9600, 9700 and 2400, control evaporation and refluxing by overlaying the mix with 50 to 100 μ L of mineral oil (Sigma Chemical Co., St. Louis, MO. Product No. M5904). The oil should not interfere when withdrawing samples. If the entire volume is to be recovered, 50 μ L of high purity chloroform should be added after amplification. The aqueous phase containing the DNA will then float on the chloroform-oil mixture, allowing easy collection. The use of AmpliWax[®] PCR Gem 50 (**Part No. N808-0150**; or **N808-0100** for 100 μ L reaction volumes) as an alternative to oil may also ease sampling and increase reproducibility in amplification when used for hot starts at low target number or with high background DNA concentrations.¹⁵

3.3 The reagent concentration ranges in the Reaction Mix recipe in **Section 2.3** are a useful starting place for amplification of different DNA targets using primers designed by the user. Optimization of reactions for each primer-template pair may be necessary and can be achieved by varying magnesium chloride concentration, primer concentration, dNTP concentration, and cycle conditions. The effect of these variations can be monitored by examining the intensity and distribution of bands after electrophoresis on agarose followed by visualization with ethidium bromide staining of the gel.

3.4 The DNA segment to be amplified from the template can be up to 10 kb long,¹⁶ although 100 to 1000 bases are more typical and easier to amplify. Start with enough copies of the template to be sure of obtaining a signal after 25 or 30 cycles; preferably $> 10^4$ copies but less than 500 ng total sample DNA per 50 μ L. Low concentrations of target DNA may require up to 35 or more cycles to produce sufficient product for analysis.

3.5 If proteases are present in the sample DNA (such as impure genomic DNA), inactivate the proteases by heating samples to 95°C for 5 minutes before adding AmpliTaq DNA Polymerase. The step can be carried out automatically with any of the GeneAmp PCR Instrument Systems.

3.6 The single-strand DNA primers should typically be 15 to 30 bases in length. The % G+C of primers should be near 50%, to maximize specificity. To avoid potential problems, primers should be purified by gel electrophoresis or HPLC ion-exchange chromatography. The optimal primer concentrations need to be determined empirically, by testing concentrations in the range of 0.1 to 1 μ M. Primer concentrations that are too low will result in little or no PCR product, while concentrations that are too high may result in amplification of non-target sequences. Primer concentrations in the range of 0.2 to 0.5 μ M should work for most PCR amplifications. Primer sequences should not complement within themselves or to each other, particularly at the 3' ends.

3.7 The optimal magnesium chloride concentration needs to be determined empirically, by testing concentrations up to 4 mM MgCl₂ for each primer set. Too little or too much MgCl₂ could reduce amplification efficiency or result in non-specific products. If the samples contain EDTA or other chelators, raise the MgCl₂ concentration in the reaction mix proportionately. Magnesium chloride concentrations should also be adjusted in parallel with significant changes in the concentration (higher or lower) of sample DNA and dNTPs.

3.8 Keep concentrations of dNTPs in the reaction mix balanced: if the concentration of any one dNTP is significantly different from the rest, AmpliTaq DNA Polymerase may tend to misincorporate them, slow down, and terminate prematurely.³

3.9 A range of volumes, dependent on the GeneAmp PCR Instrument System, will permit good amplification (see **Table 1** in **Section 4.1**). If variable volumes of reagents or template are used in the Master Mix, adjust the volume of water in the Master Mix by an equivalent amount to keep the concentrations of other reactants constant.

3.10 AmpliTaq DNA Polymerase gives optimal results under reagent conditions provided by GeneAmp 10X PCR Buffer and GeneAmp 10X Buffer II. If other PCR buffers are used, or sample preparation reagents are added that change the pH or ionic strength of the PCR, then pH and ionic strength must be optimized.

4.0 TEMPERATURE CYCLING AND CYCLE OPTIMIZATION

4.1 AmpliTaq DNA Polymerase is optimized and quality control tested for performance on GeneAmp PCR Instrument Systems. Optimal performance of the PCR process is influenced by choice of temperature, time at temperature, and length of time between temperatures for each step in the cycle.

A typical cycle consists of a melting step (94°-96°C), to separate the complementary strands of DNA; a primer annealing step (37°-65°C), depending on the target, to allow hybridization of the primers to the ssDNA and initiation of polymerization; and a primer extension step (72°C), to complete the copy initiated during annealing. Begin with a template melting step (at least 1 minute at 94°-96°C), then start cycling. Table 1 shows typical profile times for various GeneAmp PCR Instrument Systems. See their instruction manuals for details.

Samples can be stored at or below 4°C (with the GeneAmp Carryover Prevention Kit [Part No. N808-0068] use 72°C). The following steps 4.2 to 4.4 discuss the effect of changing various conditions.

Table 1. Applied Biosystems GeneAmp PCR Instrument Systems' Cycling Profile Times and Temperatures for AmpliTaq DNA Polymerase.

GeneAmp PCR Instrument System	Tubes and Volumes		Examples of Times and Temperatures			
	Tube Type	Volume Range in μ L/Tube	Initial Step	Each of 25 or more Cycles		Final Step**
				Melt	Anneal/Extend	
DNA Thermal Cycler 480	GeneAmp Thin-walled Reaction	10-100	STEP CYCLE	STEP CYCLE		TIME DELAY
			1 min 45 sec 94°-96°C 1 cycle	45 sec 94°-96°C	45 sec-1 min 37°-72°C	6 min 72°C
GeneAmp PCR System 9700* System 9600 System 2400	MicroAmp Reaction	5-100	HOLD	CYCL (CYCLE)		HOLD
			1 min 94°-96°C 1 cycle	15 sec 94°-96°C	15-30 sec 37°-72°C	6 min** 72°C

* GeneAmp PCR System 9700 depends on mode.

** After the final step, store amplified samples in refrigerator or freezer until ready to use.

4.2 High G+C content DNA needs very high annealing (> 65°C) and melting temperatures or the use of 7-deaza-2'-deoxy-GTP mixed with dGTP, to overcome secondary structures.¹⁷ The half life of AmpliTaq DNA Polymerase (< 35 minutes at > 95°C)¹⁸ suggests 95°-96°C as the maximum practical melting temperature. It is very important in the early cycles to be sure to melt completely the template DNA. When genomic DNA is used as the starting template, melting at 97°C for the first few cycles will ensure single-stranded template for the PCR reaction. The melting temperature can be reduced for the later cycles because the smaller PCR product usually melts completely at a lower temperature (unless the PCR product is excessively G+C rich) than the starting genomic DNA.

4.3 Higher annealing temperatures (> 45°C) generally result in much more specific product.^{2,19} The optimum annealing temperature can be determined empirically by testing at 5°C or smaller increments until the maximum in specificity is reached. At these temperatures AmpliTaq DNA Polymerase has significant activity so extension of primed templates is occurring. Two temperature rather than three temperature cycles may work where the anneal temperature is above 50°C; the extension is completed at the anneal temperature. The slower extension rate at lower temperatures may require optimization of the time for such combined anneal-extend steps.

4.4 The length of the target sequence will affect the required extension time. Typically, AmpliTaq DNA Polymerase has an extension rate of 2,000 to 4,000 bases per minute at 70°-80°C. Polymerization rates are significant even below 55°C and with some templates, up to 85°C.³ As the amount of DNA increases in later cycles, the number of AmpliTaq DNA Polymerase molecules may become limiting for the extension time allotted. Increasing the extend times in later cycles may be needed to maintain efficiency of amplification. Use the autosegment extension feature for the DNA Thermal Cycler 480, or the AUTO program for the GeneAmp PCR System 9600, and AUTO X for the GeneAmp PCR System 9700/2400.

TROUBLESHOOTING



<u>Observation</u>	<u>Possible Cause</u>	<u>Recommended Action</u>
Reduce or no product band is observed.	Experimental sample concentration too low.	Increase experimental sample concentration.
	Experimental sample DNA damaged or degraded.	Use experimental sample that has been processed to minimize nicking and shearing.
	Enzyme concentration too low.	Increase enzyme concentration in increments of 0.25 Unit per 50 μ L reaction.
	Magnesium concentration too low.	Increase magnesium concentration in increments of 0.2 mM.
	Denaturation time too long or too short.	Adjust denaturation time in increments of 5 seconds.
	Denaturation temperature too high or too low.	Adjust denaturation temperature in increments of 1 degree.
	Anneal/extend temperature too high.	Lower temperature in increments of 2°C.
	Anneal/extend time too short.	Lengthen time in increments of 30 seconds.
	Cycle number too low.	Increase cycle number in increments of 3.
pH or ionic strength not optimal.	See Section 3.10 .	
Primer design not optimal.	Review primer design and composition.	

<u>Observation</u>	<u>Possible Cause</u>	<u>Recommended Action</u>
Product is multi-banded.	Anneal/extend temperature too low.	Raise temperature in increments of 2 degrees.
	Primer design not optimal.	Review primer design and composition.
	Reaction Not Hot Started.	Try Manual Hot Start or AmpliTaq Gold® DNA Polymerase.
Product band is smeared.	Carry-over contamination.	Set up PCR reactions in a separate area.
	Enzyme concentration too high.	Decrease enzyme concentration in increments of 0.25 Unit per 50 µL reaction.
	Magnesium concentration too high.	Decrease magnesium concentration in increments of 0.1 mM.
	Denaturation time too short.	Lengthen denaturation time in increments of 5 seconds.
	Denaturation temperature too low.	Raise denaturation temperature in increments of 1 degree.
	Anneal/extend time too long.	Shorten time in increments of 30 seconds.
	Cycle number too high.	Reduce cycle number in increments of 3 cycles.
Partial or complete loss of product band.	Reaction Not Hot Started.	Try Manual Hot Start or AmpliTaq Gold DNA Polymerase.
	Carry-over contamination.	Set up PCR reactions in a separate area.
	Experimental sample DNA degraded.	Test a new aliquot of sample.


REFERENCES

1. Lawyer, F.C., Stoffel, S., Saiki, R.K. et al. 1989. Isolation, characterization, and expression in *E. coli* of the DNA polymerase gene from the extreme thermophile, *Thermus aquaticus*. *Journal of Biological Chemistry* **264**: 6427-6437.
2. Saiki, R.K., Gelfand, D.H., Stoffel, S. et al. 1988. Primer-directed enzymatic amplification of DNA with a thermostable DNA polymerase. *Science* **239**: 487-491.
3. Innis, M.A., Myambo, K.B., Gelfand, D.H. et al. 1988. DNA sequencing with *Thermus aquaticus* DNA polymerase and direct sequencing of polymerase chain reaction-amplified DNA. *Proceedings of National Academy of Sciences, USA*. **85**: 9436-9440.
4. Gyllensten, U. and Erlich, H.A. 1988. Generation of single stranded DNA by the polymerase chain reaction and its application to direct sequencing of the HLA DQ alpha locus. *Proceedings of the National Academy Sciences, USA*. **85**: 7652-7656.
5. Innis, M.A., Gelfand, D.H., Sninsky, J.J. et al. 1990. *PCR Protocols. A Guide to Methods and Applications*. Academic Press, Inc., San Diego, CA.
6. Kawasaki, E.S., Clark, S.C., Coyne, M.Y. et al. 1988. Diagnosis of chronic myeloid and acute lymphocytic leukemias by detection of leukemia-specific mRNA sequences amplified *in vitro*. *Proceedings of the National Academy of Sciences, USA* **85**: 5698-5702.
7. Mullis, K.B. and Faloona, F.A. 1987. Specific synthesis of DNA *in vitro* via a polymerase-catalyzed chain reaction. *Methods in Enzymology* **155**: 335-350.
8. Saiki, R.K., Scharf, S., Faloona, F. et al. 1985. Enzymatic amplification of β -globin genomic sequences and restriction site analysis for diagnosis of sickle cell anemia. *Science* **230**: 1350-1354.
9. Scharf, S.J., Horn, G.T., and Erlich, H.A. 1986. Direct cloning and sequence analysis of enzymatically amplified genomic sequences. *Science* **233**: 1076-1078.
10. Richardson, C.C. 1966. DNA polymerase from *Escherichia coli*. In: *Procedures in Nucleic Acid Research*, eds. G.L. Cantoni and D.R. Davies. Harper & Row, Inc., New York, NY, pp.263-276.
11. Gelfand, D.H., Roche Molecular Systems, Inc. 1992. Personal communication.
12. Holland, P.M., Abramson, R.D., Watson, R. et al. 1991. Detection of specific polymerase chain reaction product by utilizing the 5' -3' exonuclease activity of *Thermus aquaticus* DNA polymerase. *Proceedings of the National Academy of Sciences, USA* **88**: 7276-7280.
13. Higuchi, R. and Kwok, S. 1989. Avoiding false positives with PCR. *Nature* **339**: 237-238.

14. Kwok, S. 1990. Procedures to minimize PCR-product carry-over. In: *PCR Protocols. A Guide to Methods and Applications*, eds. Innis, M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J. Academic Press, Inc., San Diego, CA . pp. 142-145.
15. Orrego, C. 1990. Organizing a laboratory for PCR work. In: *PCR Protocols. A Guide to Methods and Applications*. eds. Innis, M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J. Academic Press, Inc., San Diego, CA. pp.447-454.
16. Chou, Q., Russell, M., Birch, D.E. et al. 1992. Prevention of pre-PCR mis-priming and primer dimerization improves low-copy-number amplifications. *Nucleic Acids Research* **20**: 1717-1723.
17. Jeffreys, A.J., Wilson, V., Neuman, R. et al. 1988. Amplification of human minisatellites by the polymerase chain reaction: Towards DNA fingerprinting of single cells. *Nucleic Acids Research* **16**: 10953-10971.
18. McConlogue, L., Brow, M.A.D., and Innis, M.A. 1988. Structure-independent DNA amplification by PCR using 7-deaza-2'-deoxyguanosine. *Nucleic Acids Research* **16**: 9869.
19. Gelfand, D.H., and White, T.J. 1990. Thermostable DNA Polymerases. In: *PCR Protocols. A Guide to Methods and Applications*. eds. Innis, M.A., Gelfand, D.H., Sninsky, J.J., and White, T.J., Academic Press, Inc., San Diego, CA, pp.129-141.
20. Rychlik, W., Spencer, W.J., and Rhoads, R.E. 1990. Optimization of the annealing temperature for DNA Amplification *in vitro*. *Nucleic Acids Research* **18**: 6409-6412.



Applied Biosystems



ORDERING INFORMATION
PRODUCTS ARE NOT FOR USE IN DIAGNOSTIC PROCEDURES

In the United States and Canada:	Outside of the United States:
For PCR Technical Support call toll-free 1 800 762 4001 FAX: 240 453 4613	For the number of your local Applied Biosystems Sales Office call:
To order PCR consumables and reagents call toll-free 1 800 327 3002	Australia 61 3 9730 8600
To order PCR instruments or to reach your sales representative call toll-free 1 800 345 5224	Austria 43 1 867 35 75 0
For PCR product literature call toll-free 1 800 345 5224	Belgium 32 (0) 2 712 555
	Czech Rep. & Slovakia 420 2 61 222 164
	Denmark 45 45 58 6000
	Finland 358 (0) 9 251 24 250
	France 33 (0) 1 69 59 85 85
	Germany 49 (0) 6150 10 10
	Hungary 36 (0) 1 270 8398
	Italy 39 039 83891
	Japan 81-3-5566-6100
	Latin America - Mexico 52 5 651 7077
	Middle East/North Africa 39 (0) 39 8389 481
	Netherlands (Benelux) 31 (0) 180 331 400
	Norway 47 23 12 06 05
	Poland 48 22 866 40 10
	Portugal 351 (0) 22 605 33 14
	South Africa 27 11 4780411
	South East Europe 385 1 34 91 927
	Spain 34 91 806 1210
	Sweden 46 (0) 8 619 4400
	Switzerland 41 (0) 41 799 7777
	United Kingdom 44 (0) 1925 825650
	All other countries not listed 44 (0) 1925 282481

GeneAmp® PCR Instrument Systems from Applied Biosystems:

N801-0100 DNA Thermal Cycler 480, 120 V, 60 Hz

N801-0101 DNA Thermal Cycler 480, 240 V, 50 Hz

N801-0102 DNA Thermal Cycler 480, 100 V, 50/60 Hz

N803-0001 GeneAmp® PCR System 2400, 120 V, 60 Hz

N803-0002 GeneAmp® PCR System 2400, 240 V, 50 Hz

N803-0003 GeneAmp® PCR System 2400, 100 V, 50/60 Hz

N801-0001 GeneAmp® PCR System 9600, 120 V, 60 Hz

N801-0002 GeneAmp® PCR System 9600, 240 V, 50 Hz

N801-0003 GeneAmp® PCR System 9600, 100 V, 50/60 Hz

N805-0001 GeneAmp® PCR System 9700, Multi-voltage

AmpliCycle™ is a trademark and AmpliTaq®, AmpliTaq Gold®, GeneAmp® and AmpliWax® are registered trademarks of Roche Molecular Systems, Inc.

MicroAmp® is a registered trademark and Applied Biosystems is a trademark of Applied Biosystems.

Nonidet® P 40 is a registered trademark of Shell Oil Company.

Roche® is a registered trademark of Hoffmann-La Roche Inc.

Tween® is a registered trademark of ICI.

Applied Biosystems

This product is manufactured for Applied Biosystems by Roche Molecular Systems, Inc., Branchburg, New Jersey USA



Applied Biosystems
850 Lincoln Centre Drive
Foster City, CA 94404