

CsCl₂ Maxi DNA Prep

Experimental Design Considerations

- Yield depends on bacterial strain (DH5 typically gives high yields while BL21 gives poor yields)
- Yield also depends on backbone of plasmid with old pBR322 plasmids and GEX plasmids giving low yields
- Better yields are obtained if cells are allowed to grow at saturation for several hours. Therefore, cells should be cultured for at least the minimum incubation times listed

Protocol

- Inoculate 5ml media (+100ug/ml Amp) in Morning
- Inoculate 250ml 2XYT (100ug/ml Amp) with 5ml culture in afternoon.
- Incubate O/N (at least 16hrs).
- Spin 5 minutes at 5K in 500ml centrifuge bottles
- Suspend pellet in 6ml Soln. I and transfer to 50ml conical tube (polypropylene).
- Add 12 ml Soln. II (0.2N NaOH, 1% SDS), shake vigorously.
- Add 9 ml Soln. III, rock back and forth GENTLY until you see white fluffy precipitate.
- Spin 5 min. on high (ca. 3K) in clinical centrifuge.
- Transfer aqueous to new 50 ml tube (don't transfer white wafer on top).
- Add 10 ml Phenol/CHCl₃, and shake vigorously
- Spin at high speed for 5 minutes
- Transfer top phase to new 50 ml tube.
- Add 100% EtOH up to 50ml.
- Let sit ca. 20-40 min.
- Spin 20 min. clinical centrifuge (high speed).
- Pour off EtOH, wipe out excess EtOH with kimwipe.
- Add 4ml TE, disrupt pellet so that it floats on surface - it will then dissolve on it's own.
- Add 4.7g CsCl₂.
- Add 200ul 10mg/ml Ethidium Bromide.
- Transfer to vTI80 opti-seal centrifuge tubes.
- Spin 4hrs, 75k or O/N at 65k.
- Take out DNA with syringe (21-gauge needle), take off needle before transferring to 15ml polystyrene conical tube.
- Extract several times with 1ml of H₂O saturated n-butanol until clear [if there is difficulty extracting EtBr, it probably means that the EtBr is old].
- Add TE to bring up volume to 3ml.
- Add 300ul 3M NaOAc, +8ml EtOH, let sit at RT 20 min (if you have low yield, keep at -20C overnight)
- Spin in clinical centrifuge 20min. on high.
- Wash with cold 70% EtOH.
- Let dry [we usually wrap two kimwipes around the eraser end of a pencil to dry the inside of tube (2X)]
- Resuspend in TE (typically 500-1000ul)

Buffers

SOLUTION I

50mM glucose	<u>1L</u> 9g
25mM Tris (8.0)	25ml 1M
10mM EDTA	20ml 0.5M

SOLUTION III

5M KOAc
Glacial Acetic Acid

1 liter

295g KOAc
115ml Acetic Acid

2xYT

tryptone	<u>1L</u> 16g
yeast extract	10g
NaCl	5g

TE

10mM Tris (pH 7.6)
1mM EDTA (pH 8.0)