

## Product Contents

### GoTaq® Hot Start Polymerase:

#### Supplied With:

Cat. #	GoTaq® Hot Start Polymerase	5X Green GoTaq® Flexi Buffer	5X Colorless GoTaq® Flexi Buffer	Magnesium Chloride Solution, 25mM
M5001	100 units (M500A)	1ml (M891A)	1ml (M890A)	0.75ml (A351B)
M5005	500 units (M500B)	4 × 1ml (M891A)	4 × 1ml (M890A)	3 × 1.2ml (A351H)

**Description:** GoTaq® Hot Start Polymerase(a-c) contains GoTaq® Hot Start Polymerase, 5X Green GoTaq® Flexi Buffer, 5X Colorless GoTaq® Flexi Buffer and 25mM MgCl<sub>2</sub>. GoTaq® Hot Start Polymerase contains the high-performance GoTaq® DNA polymerase bound to a proprietary antibody that blocks polymerase activity. The polymerase activity is restored during the initial denaturation step when amplification reactions are heated at 94–95°C for two minutes. This allows hot-start PCR in which polymerase activity is inhibited at temperatures below 70°C, allowing convenient, room-temperature reaction setup. Hot-start PCR is advantageous for some amplification targets, because it may eliminate or minimize primer-dimer and secondary products. In some cases, hot-start PCR may improve yields. GoTaq® Hot Start Polymerase exhibits 5'→3' exonuclease activity.

The 5X Green GoTaq® Flexi Buffer contains two dyes (a blue dye and a yellow dye) that separate during electrophoresis to monitor migration progress. The colorless buffer is used when direct fluorescence or absorbance readings are required without prior purification of the amplified DNA from the PCR. The Flexi Buffers do not contain magnesium, allowing easy optimization in amplification reactions.

**Biological Source:** The enzyme is derived from bacteria. The antibody is derived from murine cell culture.

**Enzyme Concentration:** 5u/μl.

**5X Green GoTaq® Flexi Buffer (Part # M891A):** Proprietary formulation supplied at pH 8.5 containing blue dye and yellow dye. The blue dye migrates at the same rate as a 3–5kb DNA fragment in a 1% agarose gel. The yellow dye migrates at a rate faster than primers (<50bp) in a 1% agarose gel. Green GoTaq® Flexi Buffer also increases the density of the sample, so it will sink into the well of the agarose gel, allowing reactions to be loaded directly onto gels without loading dye. This buffer does not contain magnesium.

**5X Colorless GoTaq® Flexi Buffer (Part # M890A):** Proprietary formulation supplied at pH 8.5. This buffer does not contain magnesium.

**Magnesium Chloride Solution, 25mM (Part# A351B, A351H):** Provided to allow users to optimize MgCl<sub>2</sub> concentration according to their individual requirements. **Vortex the MgCl<sub>2</sub> thoroughly after thawing and prior to use.**

**Storage Conditions:** See the Product Information Label for storage conditions. Avoid exposure to frequent temperature changes. See the expiration date on the Product Information label.

**Unit Definition:** One unit is defined as the amount of enzyme required to catalyze the incorporation of 10 nanomoles of dNTPs into acid-insoluble material in 30 minutes at 74°C.

## Quality Control Assays

**Functional Assay:** GoTaq® Hot Start Polymerase is tested for performance in PCR using 1.25 units of enzyme to amplify a 360bp region of the α-1-antitrypsin gene and a 2.4kb region of the APC gene from 100 molecules (0.35ng) of human genomic DNA in separate reactions. The resulting PCR products are visualized as single bands on an ethidium bromide-stained agarose gel.

**Nuclease Assays:** No contaminating endonuclease or exonuclease activity detected.

**Hot-Start Amplification Assay:** GoTaq® Hot Start Polymerase is tested in PCR for its ability to amplify a hot-start model template and produce a single 1.5kb band, eliminating extraneous bands. In a non-hot-start PCR, this template produces an additional band at 410bp.

**5X Green GoTaq® Flexi Buffer Migration Pattern:** The 5X Green GoTaq® Flexi Buffer does not interfere with the migration of a 1kb DNA ladder when it is used as a loading dye for agarose gel electrophoresis.



#### PCR Satisfaction Guarantee

Promega's PCR Systems, enzymes and reagents are proven in PCR to ensure reliable, high performance results. Your success is important to us. Our products are backed by a worldwide team of Technical Support scientists. Please contact them for applications or technical assistance. If you are not completely satisfied with any Promega PCR product we will send a replacement or refund your account.  
*That's Our PCR Guarantee!*

Product must be within expiration date and have been stored and used in accordance with product literature. See Promega Product Insert for specific tests performed.

(a)Use of this product for basic PCR is outside of any valid US or European patents assigned to Hoffman La-Roche or Applera. This product can be used for basic PCR in research without any license or royalty fees. This product is for research use only.

(b)U.S. Pat. No. 6,242,235, Australian Pat. No. 761757, Canadian Pat. No. 2,335,153, Chinese Pat. No. ZL99808861.7, Hong Kong Pat. No. HK 1040262, Japanese Pat. No. 3673175 and European Pat. No. 1088060 have been issued to Promega Corporation for enzyme stabilization by cationic surfactants. Other patents are pending.

(c)Licensed under U.S. Pat. Nos. 5,338,671 and 5,587,287 and corresponding patents in other countries.

Part# 9PIM500

Revised 7/08



## Promega

#### Promega Corporation

2800 Woods Hollow Road	
Madison, WI 53711-5399	USA
Telephone	608-274-4330
Toll Free	800-356-9526
Fax	608-277-2516
Internet	www.promega.com

#### PRODUCT USE LIMITATIONS, WARRANTY DISCLAIMER

Promega manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products. Promega products contain chemicals which may be harmful if misused. Due care should be exercised with all Promega products to prevent direct human contact.

Each Promega product is shipped with documentation stating specifications and other technical information. Promega products are warranted to meet or exceed the stated specifications. Promega's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Promega makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO PROMEGA PRODUCTS. In no event shall Promega be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of Promega products to perform in accordance with the stated specifications.

© 2008 Promega Corporation. All Rights Reserved.

GoTaq and Wizard are registered trademarks of Promega Corporation.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.

Part# 9PIM500

Printed in USA. Revised 7/08

## I. Standard Application

### Reagents to Be Supplied by the User

- PCR Nucleotide Mix (Cat.# C1141)
  - upstream primer
  - template DNA
  - Nuclease-Free Water (Cat.# P1193)
  - downstream primer
  - mineral oil (optional)
- In a sterile, nuclease-free microcentrifuge tube, combine the following components at room temperature:

Component	Final Volume	Final Conc.
5X Green or Colorless GoTaq® Flexi Buffer <sup>1</sup>	10µl	1X
MgCl <sub>2</sub> Solution, 25mM <sup>1</sup>	2–8µl	1.0–4.0mM
PCR Nucleotide Mix, 10mM each	1µl	0.2mM each dNTP
upstream primer	Xµl	0.1–1.0µM
downstream primer	Yµl	0.1–1.0µM
GoTaq® Hot Start Polymerase (5u/µl)	0.25µl	1.25u
template DNA	Zµl	<0.5µg/50µl
Nuclease-Free Water to	50µl	

<sup>1</sup>Thaw completely and vortex thoroughly prior to use.

- If using a thermal cycler **without a heated lid**, overlay the reaction mix with 1–2 drops (approximately 50µl) of mineral oil to prevent evaporation during thermal cycling. Centrifuge the reactions in a microcentrifuge for 5 seconds.
- Reactions can be placed in a room-temperature thermal cycler. **A 2-minute initial denaturation step at 94–95°C is required to inactivate the antibody and initiate hot-start PCR.**

## II. General Guidelines for Amplification by PCR

### A. Denaturation

- Following the initial 2-minute 94–95°C denaturation step, denaturation steps should be between 30 seconds and 1 minute.

### B. Annealing

- Optimize the annealing conditions by performing the reaction starting approximately 5°C below the calculated melting temperature of the primers and increasing the temperature in increments of 1°C to the annealing temperature.
- The annealing step is typically 30 seconds to 1 minute.

### C. Extension

- The extension reaction is typically performed at the optimal temperature for Taq DNA polymerase, which is 72–74°C.
- Allow approximately 1 minute for every 1kb of DNA to be amplified.
- A final extension of 5 minutes at 72–74°C is recommended.

### D. Refrigeration

- If the thermal cycler has a refrigeration or “soak” cycle, the cycling reaction can be programmed to end by holding the tubes at 4°C for several hours.
- This cycle can minimize any polymerase activity that might occur at higher temperatures, although this is not usually a problem.

### E. Cycle Number

- Generally, 25–30 cycles result in optimal amplification of desired products.
- Occasionally, up to 40 cycles may be performed, especially for detection of low-copy targets.

## II. General Considerations

### A. Buffer Choice

We recommend using the 5X Green GoTaq® Flexi Buffer in any amplification reaction that will be visualized by agarose gel electrophoresis followed by ethidium-bromide staining. The 5X Green GoTaq® Flexi Buffer is not recommended for any downstream applications using absorbance or fluorescence excitation, because the yellow and blue dyes in the

reaction buffer may interfere with these applications. The dyes absorb between 225–300nm, making standard A<sub>260</sub> readings to determine DNA concentration unreliable. Also, the dyes have excitation peaks at 488nm and between 600–700nm that correspond to the excitation wavelengths commonly used in fluorescence detection instrumentation. However, for some instrumentation, such as a fluorescent gel scanner that uses a 488nm excitation wavelength, there will be minimal interference, since it is the yellow dye that absorbs this wavelength. Gels scanned by this method will have a light gray dye front below the primers corresponding to the yellow dye front. The Green and Colorless GoTaq® Flexi Buffers give approximately equivalent amplification yields. Balanced amplifications between the two buffers may require further optimization.

For reactions going directly from thermal cycler to an application using absorbance or fluorescence, the 5X Colorless GoTaq® Flexi Buffer is recommended. If both agarose gel analysis and further downstream applications involving absorbance or fluorescence will be used, the two dyes can be removed from the Green GoTaq® Flexi reactions using standard PCR clean-up systems like the Wizard® SV Gel and PCR Clean-Up System (Cat.# A9281) or the Wizard® SV 96 PCR Clean-Up System (Cat.# A9341).

### B. Enzyme Concentration

Promega has found that 1.25 units of GoTaq® DNA Polymerase per 50µl amplification reaction is adequate for most amplifications. Adding extra enzyme generally does not produce significant increases in yield. However, in some cases, more enzyme may be beneficial.

### C. Primer Design

PCR primers generally range in length from 15–30 bases and are designed to flank the region of interest. Primers should contain 40–60% (G + C), and care should be taken to avoid sequences that might produce internal secondary structure. The 3'-ends of the primers should not be complementary to avoid the production of primer-dimers. Primer-dimers unnecessarily deplete primers from the reaction and result in an unwanted polymerase reaction that competes with the desired reaction. Avoid three G or C nucleotides in a row near the 3'-end of the primer, because this may result in nonspecific primer annealing, increasing the synthesis of undesirable reaction products. Ideally, both primers should have nearly identical melting temperatures (T<sub>m</sub>); in this manner, the two primers should anneal roughly at the same temperature. The annealing temperature of the reaction depends on the primer with the lowest T<sub>m</sub>. For assistance with calculating the T<sub>m</sub> of any primer, a T<sub>m</sub> Calculator is provided on the BioMath page of the Promega web site at: [www.promega.com/biomath/](http://www.promega.com/biomath/)

### D. Amplification Troubleshooting

To overcome low yield or no yield in amplifications, we recommend the following suggestions:

- Adjust annealing temperature. The reaction buffer composition affects the melting properties of DNA. See BioMath Calculator to calculate the melting temperature for primers in the GoTaq® reaction ([www.promega.com/biomath/](http://www.promega.com/biomath/)).
- Minimize the effect of amplification inhibitors. Some DNA isolation procedures, particularly genomic DNA isolation, can result in the copurification of amplification inhibitors. Reduce the volume of template DNA in reaction, or dilute template DNA prior to adding to reaction. Diluting samples up to 1:10,000 has been shown to be effective in improving results, depending on initial DNA concentration.
- Increase template DNA purity. Include an ethanol precipitation and wash step prior to amplification to remove inhibitors that copurify with the DNA.
- Add PCR additives. Adding PCR-enhancing agents (e.g., DMSO or betaine) may improve yields. General stabilizing agents such as BSA (Sigma Cat.# A7030; final concentration 0.16mg/ml) also may help to overcome amplification failure.

### E. More Information on Amplification

More information on amplification is available online at the Promega web site:

PCR Applications: [www.promega.com/paguide/chap1.html](http://www.promega.com/paguide/chap1.html)

PCR Core Systems Technical Bulletin: [www.promega.com/tbs/tb254/tb254.html](http://www.promega.com/tbs/tb254/tb254.html)

DNA Analysis Notebook: [www.promega.com/guides/dna\\_guide/](http://www.promega.com/guides/dna_guide/)