

pGEM[®]-T and pGEM[®]-T Easy Vector Systems

Technical Manual No. 042

INSTRUCTIONS FOR USE OF PRODUCTS A1360, A1380, A3600 AND A3610.

PLEASE DISCARD PREVIOUS VERSIONS.

All Technical Literature is Available on the Internet at www.promega.com

Please visit the web site to verify that you are using the most current version of this Technical Manual.

I. Description	2
II. Vector Maps.....	3
A. pGEM [®] -T Vector and pGEM [®] -T Easy Vector Multiple Cloning Sequences	3
B. pGEM [®] -T Vector Map and Sequence Reference Points	4
C. pGEM [®] -T Easy Vector Map and Sequence Reference Points	5
III. Product Components	6
IV. Protocol for Ligations Using the pGEM [®] -T and pGEM [®] -T Easy Vectors and the 2X Rapid Ligation Buffer	7
V. Protocol for Transformations Using the pGEM [®] -T and pGEM [®] -T Easy Vector Ligation Reactions.....	7
VI. General Considerations	9
A. PCR Product Purity	9
B. Blunt-Ended PCR Products	9
C. Optimizing Insert:Vector Molar Ratios	11
D. Screening Transformants for Inserts.....	12
E. Experimental Controls	12
VII. Isolation of Recombinant Plasmid DNA	13
VIII. Generation of Single-Stranded DNA from the pGEM [®] -T and pGEM [®] -T Easy Vectors	14
IX. Troubleshooting.....	14
X. References	18
XI. Appendix A: Vector Sequences and Restriction Sites	18
A. pGEM [®] -T Vector Sequence.....	18
B. pGEM [®] -T Vector Restriction Sites.....	20
C. pGEM [®] -T Easy Vector Sequence.....	22
D. pGEM [®] -T Easy Vector Restriction Sites.....	24
XII. Appendix B: Reference Information.....	25
A. Composition of Buffers and Solutions	25
B. Related Products	26

I. Description

The pGEM[®]-T and pGEM[®]-T Easy Vector Systems^(a,b) are convenient systems for the cloning of PCR products. The vectors are prepared by cutting Promega's pGEM[®]-5Zf(+)^(b) and pGEM[®]-T Easy Vectors with *EcoR* V and adding a 3' terminal thymidine to both ends. These single 3'-T overhangs at the insertion site greatly improve the efficiency of ligation of a PCR product into the plasmids by preventing recircularization of the vector and providing a compatible overhang for PCR products generated by certain thermostable polymerases (1,2). As summarized in Table 1, these polymerases often add a single deoxyadenosine, in a template-independent fashion, to the 3'-ends of the amplified fragments (3,4).

The high copy number pGEM[®]-T and pGEM[®]-T Easy Vectors contain T7 and SP6 RNA polymerase promoters flanking a multiple cloning region within the α -peptide coding region of the enzyme β -galactosidase. Insertional inactivation of the α -peptide allows recombinant clones to be directly identified by color screening on indicator plates. The multiple cloning region of the two vectors includes restriction sites conveniently arranged for use with Promega's Erase-a-Base[®] System (Cat.# E5750) for generating nested sets of deletions.

Both the pGEM[®]-T and pGEM[®]-T Easy Vector contain multiple restriction sites within the multiple cloning region. These restriction sites allow for the release of the insert by digestion with a single restriction enzyme. The pGEM[®]-T Easy Vector multiple cloning region is flanked by recognition sites for the restriction enzymes *EcoR* I, *BstZ* I and *Not* I, thus providing three single-enzyme digestions for release of the insert, while the pGEM[®]-T Vector cloning region is flanked by recognition sites for the enzyme *BstZ* I. Alternatively, a double-digestion may be used to release the insert from either vector.

The pGEM[®]-T and pGEM[®]-T Easy Vectors also contain the origin of replication of the filamentous phage f1 for the preparation of single-stranded DNA (ssDNA; see Section VII). The ssDNA molecule exported corresponds to the bottom strand shown in Figure 1, Panels A and B, for the pGEM[®]-T and pGEM[®]-T Easy Vectors (nonrecombinant), respectively.

The pGEM[®]-T and pGEM[®]-T Easy Vector Systems include a 2X Rapid Ligation Buffer for ligation of PCR products. Reactions using this buffer may be incubated for 1 hour at room temperature. The incubation period may be extended to increase the number of colonies after transformation. Generally, an overnight incubation at 4°C will produce the maximum number of transformants.


Table 1. Comparison of PCR Product Properties for Some Thermostable DNA Polymerases.

Characteristic	Thermostable DNA Polymerase						
	<i>Taq</i> AmpliTaq [®]	<i>Tfl</i>	<i>Tth</i>	<i>Vent</i> [®] / (<i>Tli</i>)	<i>Deep</i> <i>Vent</i> [®]	<i>Pfu</i>	<i>Pwo</i>
Resulting DNA ends	3' A	3' A	3' A	>95% Blunt	>95% Blunt	Blunt	Blunt
5'→3' exonuclease activity	Yes	Yes	Yes	No	No	No	No
3'→5' exonuclease activity	No	No	No	Yes	Yes	Yes	Yes

II. Vector Maps


A. pGEM[®]-T Vector and pGEM[®]-T Easy Vector Multiple Cloning Sequences

pGEM[®]-T Vector



03571MA06/2A

pGEM[®]-T Easy Vector



1517MA06/6A

Note: A single digest with *BstZ I* (Cat.# R6881) will release inserts cloned into the pGEM[®]-T Vector.

Note: A single digest with *BstZ I* (Cat.# R6881), *EcoR I* (Cat.# R6011) or *Not I* (Cat.# R6431) will release inserts cloned into the pGEM[®]-T Easy Vector.

Figure 1. The promoter and multiple cloning sequence of the pGEM[®]-T and pGEM[®]-T Easy Vectors. The top strand of the sequence shown corresponds to the RNA synthesized by T7 RNA polymerase. The bottom strand corresponds to the RNA synthesized by SP6 RNA polymerase.

B. pGEM[®]-T Vector Map and Sequence Reference Points

Note: Inserts can be sequenced using the following primers:
 SP6 Promoter Primer (Cat.# Q5011),
 T7 Promoter Primer (Cat.# Q5021),
 pUC/M13 Forward Primer (Cat.# Q5601),
 pUC/M13 Reverse Primer (Cat.# Q5421).




Figure 2. pGEM[®]-T Vector circle map and sequence reference points.

Note: A single digest with *BstZ* I (Cat.# R6881) will release inserts cloned into the pGEM[®]-T Vector. Double digests can also be used to release inserts.

pGEM[®]-T Vector sequence reference points:

T7 RNA polymerase transcription initiation site	1
multiple cloning region	10–113
SP6 RNA polymerase promoter (–17 to +3)	124–143
SP6 RNA polymerase transcription initiation site	126
pUC/M13 Reverse Sequencing Primer binding site	161–177
<i>lacZ</i> start codon	165
<i>lac</i> operator	185–201
β-lactamase coding region	1322–2182
phage f1 region	2365–2820
<i>lac</i> operon sequences	2821–2981, 151–380
pUC/M13 Forward Sequencing Primer binding site	2941–2957
T7 RNA polymerase promoter (–17 to +3)	2984–3

03561/A/04/3A

C. pGEM[®]-T Easy Vector Map and Sequence Reference Points




Figure 3. pGEM[®]-T Easy Vector circle map and sequence reference points.

pGEM[®]-T Easy Vector sequence reference points:

T7 RNA polymerase transcription initiation site	1
multiple cloning region	10–128
SP6 RNA polymerase promoter (–17 to +3)	139–158
SP6 RNA polymerase transcription initiation site	141
pUC/M13 Reverse Sequencing Primer binding site	176–197
<i>lacZ</i> start codon	180
<i>lac</i> operator	200–216
β -lactamase coding region	1337–2197
phage f1 region	2380–2835
<i>lac</i> operon sequences	2836–2996, 166–395
pUC/M13 Forward Sequencing Primer binding site	2949–2972
T7 RNA polymerase promoter (–17 to +3)	2999–3

Specialized applications of the pGEM[®]-T and pGEM[®]-T Easy Vectors:

- Cloning PCR products.
- Construction of unidirectional nested deletions with the Erase-a-Base[®] System.
- Production of ssDNA.
- Blue/white screening for recombinants.
- In vitro transcription from dual opposed promoters. (For protocol information, please request Promega's *Riboprobe[®] in vitro Transcription Systems^(c) Technical Manual #TM016.*)

Note: Inserts can be sequenced using the following primers:
 SP6 Promoter Primer (Cat.# Q5011),
 T7 Promoter Primer (Cat.# Q5021),
 pUC/M13 Forward Primer (Cat.# Q5601),
 pUC/M13 Reverse Primer (Cat.# Q5421).

Note: A single digest with *BstZ I* (Cat.# R6881) will release inserts cloned into the pGEM[®]-T Vector. Double digests can also be used to release inserts.

1473VA06/6A

III. Product Components

Product	Size	Cat.#
pGEM®-T Vector System I	20 reactions	A3600

For Laboratory Use. Includes:

- 1.2µg pGEM®-T Vector (50ng/µl)
- 12µl Control Insert DNA (4ng/µl)
- 100u T4 DNA Ligase
- 200µl 2X Rapid Ligation Buffer, T4 DNA Ligase
- 1 Protocol

Product	Size	Cat.#
pGEM®-T Vector System II	20 reactions	A3610

For Laboratory Use. Includes:

- 1.2µg pGEM®-T Vector (50ng/µl)
- 12µl Control Insert DNA (4ng/µl)
- 100u T4 DNA Ligase
- 200µl 2X Rapid Ligation Buffer, T4 DNA Ligase
- 1.2ml JM109 Competent Cells, High Efficiency (6 × 200µl)
- 1 Protocol

Product	Size	Cat.#
pGEM®-T Easy Vector System I	20 reactions	A1360

For Laboratory Use. Includes:

- 1.2µg pGEM®-T Easy Vector (50ng/µl)
- 12µl Control Insert DNA (4ng/µl)
- 100u T4 DNA Ligase
- 200µl 2X Rapid Ligation Buffer, T4 DNA Ligase
- 1 Protocol

Product	Size	Cat.#
pGEM®-T Easy Vector System II	20 reactions	A1380

For Laboratory Use. Includes:

- 1.2µg pGEM®-T Easy Vector (50ng/µl)
- 12µl Control Insert DNA (4ng/µl)
- 100u T4 DNA Ligase
- 200µl 2X Rapid Ligation Buffer, T4 DNA Ligase
- 1.2ml JM109 Competent Cells, High Efficiency (6 × 200µl)
- 1 Protocol

Storage Conditions: For Cat.# A3610, A1380, store the Competent Cells at -70°C. All other components can be stored at -20°C or -70°C.

IV. Protocol for Ligations Using the pGEM[®]-T and pGEM[®]-T Easy Vectors and the 2X Rapid Ligation Buffer

1. Briefly centrifuge the pGEM[®]-T or pGEM[®]-T Easy Vector and Control Insert DNA tubes to collect contents at the bottom of the tubes.
2. Set up ligation reactions as described below. **Note:** Use 0.5ml tubes known to have low DNA-binding capacity (e.g., VWR Cat.# 20170-310).
3. Vortex the 2X Rapid Ligation Buffer vigorously before each use.

	Standard Reaction	Positive Control	Background Control
2X Rapid Ligation Buffer, T4 DNA Ligase	5µl	5µl	5µl
pGEM [®] -T or pGEM [®] -T Easy Vector (50ng)	1µl	1µl	1µl
PCR product	Xµl*	–	–
Control Insert DNA	–	2µl	–
T4 DNA Ligase (3 Weiss units/µl)	1µl	1µl	1µl
deionized water to a final volume of	10µl	10µl	10µl

*Molar ratio of PCR product:vector may require optimization (see Section VI.C).

4. Mix the reactions by pipetting. Incubate the reactions 1 hour at room temperature.

Alternatively, if the maximum number of transformants is required, incubate the reactions overnight at 4°C.

Notes:

1. Use only Promega T4 DNA Ligase supplied with this system in performing pGEM[®]-T and pGEM[®]-T Easy Vector ligations. Other commercial preparations of T4 DNA ligase may contain exonuclease activities that may remove the terminal deoxythymidines from the vector.
2. 2X Rapid Ligation Buffer contains ATP, which degrades during temperature fluctuations. Avoid multiple freeze-thaw cycles and exposure to frequent temperature changes by making single-use aliquots of the buffer.
3. It is important to vortex the 2X Rapid Ligation Buffer before each use.
4. Longer incubation times will increase the number of transformants. Generally, incubation overnight at 4°C will produce the maximum number of transformants.

V. Protocol for Transformations Using the pGEM[®]-T and pGEM[®]-T Easy Vector Ligation Reactions

Use high-efficiency competent cells ($\geq 1 \times 10^8$ cfu/µg DNA) for transformations. The ligation of fragments with a single-base overhang can be inefficient, so it is essential to use cells with a transformation efficiency of 1×10^8 cfu/µg DNA (or higher) in order to obtain a reasonable number of colonies (see Section VI.E).



JM109 cells

should be maintained on M9 minimal medium plates, supplemented with thiamine hydrochloride, prior to the preparation of competent cells



In Step 3, avoid excessive pipetting as the competent cells are extremely fragile.

We recommend using JM109 High Efficiency Competent Cells (Cat.# L2001); these cells are provided with the pGEM®-T and pGEM®-T Easy Vector Systems II. Other host strains may be used, but they should be compatible with blue/white color screening and standard ampicillin selection.

JM109 cells should be maintained on M9 minimal medium plates supplemented with thiamine hydrochloride prior to the preparation of competent cells. This selects for the presence of the F' episome, containing both the *proAB* genes, which complement proline auxotrophy in a host with a (*proAB*) deletion, and *lacI^qΔM15*, required in blue/white color screening. If you are using competent cells other than JM109 High Efficiency Competent Cells purchased from Promega, it is important that the appropriate transformation protocol be followed. Selection for transformants should be on LB/ampicillin/IPTG/X-Gal plates (see Section XI.A). For best results, do not use plates that are more than 1 month old.

The genotype of JM109 is *recA1, endA1, gyrA96, thi, hsdR17 (r_K⁻,m_K⁺), relA1, supE44, Δ(lac-proAB), [F', traD36, proAB, lacI^qΔM15]* (5).

Materials to Be Supplied by the User

(Solution compositions are provided in Section XII.A.)

- LB plates with ampicillin/IPTG/X-Gal
 - SOC medium
1. Prepare 2 LB/ampicillin/IPTG/X-Gal plates for each ligation reaction, plus two plates for determining transformation efficiency (see Section VI.E). Equilibrate the plates to room temperature prior to plating (Step 10).
 2. Centrifuge the tubes containing the ligation reactions to collect contents at the bottom of the tube. Add 2μl of each ligation reaction to a sterile 1.5ml microcentrifuge tube on ice (see Note 1). Set up another tube on ice with 0.1ng uncut plasmid for determination of the transformation efficiency of the competent cells (see Section VI.E).
 3. Remove tube(s) of frozen JM109 High Efficiency Competent Cells from -70°C storage and place in an ice bath until just thawed (about 5 minutes). Mix the cells by **gently** flicking the tube.
 4. **Carefully** transfer 50μl of cells into each tube prepared in Step 2 (100μl cells for determination of transformation efficiency).
 5. **Gently** flick the tubes to mix and place them on ice for 20 minutes.
 6. Heat-shock the cells for 45–50 seconds in a water bath at exactly 42°C (**Do Not Shake**).
 7. Immediately return the tubes to ice for 2 minutes.
 8. Add 950μl room temperature SOC medium to the tubes containing cells transformed with ligation reactions and 900μl to the tube containing cells transformed with uncut plasmid (LB broth may be substituted, but colony number may be lower).
 9. Incubate for 1.5 hours at 37°C with shaking (~150rpm).
 10. Plate 100μl of each transformation culture onto duplicate LB/ampicillin/IPTG/X-Gal plates. For the transformation control, a 1:10 dilution with SOC medium is recommended for plating. If a higher number of colonies is desired, the cells may be pelleted by centrifugation at 1,000 × *g* for 10 minutes, resuspended in 200μl of SOC medium, and 100μl plated on each of 2 plates.

11. Incubate the plates overnight (16–24 hours) at 37°C. In our experience, approximately 100 colonies per plate are routinely seen when using competent cells that are 1×10^8 cfu/ μ g DNA, if 100 μ l is plated. Longer incubations or storage of plates at 4°C (after 37°C overnight incubation) may be used to facilitate blue color development. White colonies generally contain inserts; however, inserts may also be present in blue colonies. Please see Section VI.D for more information.

Notes:

1. In our experience, the use of larger (17 × 100mm) polypropylene tubes (e.g., Falcon Cat.# 2059) has been observed to increase transformation efficiency. Tubes from some manufacturers bind DNA, thereby decreasing the colony number, and should be avoided.
2. Colonies containing β -galactosidase activity may grow poorly relative to cells lacking this activity. After overnight growth, the blue colonies may be smaller than the white colonies, which are approximately one millimeter in diameter.

VI. General Considerations**A. PCR Product Purity**

An aliquot of the PCR reaction should be analyzed on an agarose gel before use in the ligation reaction. The PCR product to be ligated can be gel-purified or purified directly from the PCR amplification reaction using the Wizard® SV Gel and PCR Clean-Up System^(d) (Cat.# A9281). Exposure to shortwave ultraviolet light should be minimized in order to avoid the formation of pyrimidine dimers. If smearing of the PCR product or inappropriate banding is observed on the gel, excise the bands to be cloned and purify the DNA with Wizard® SV Gel and PCR Clean-Up System. Even if distinct bands of expected size are observed, primer-dimers should be removed by gel purification or by using the Wizard® SV Gel and PCR Clean-Up System to purify the bands of interest directly from the reaction mix. Use of crude PCR product may produce successful ligations in some cases; however, the number of white colonies containing the relevant insert may be reduced due to preferential incorporation of primer-dimers or other extraneous reaction products. Therefore, it may be necessary to screen numerous colonies in order to identify clones that contain the PCR product of interest.

B. Blunt-Ended PCR Products

Thermostable DNA polymerases with proofreading activity, such as *Pfu* DNA Polymerase^(e) (Cat.# M7741), *Pwo* DNA polymerase and *Tli* DNA Polymerase^(e) (Cat.# M7101) generate blunt-ended fragments during PCR amplification. Nevertheless, PCR fragments generated using these polymerases can be modified using the A-tailing procedure (Figure 4) and ligated into the pGEM®-T and pGEM®-T Easy Vectors (6). Using this method, only one insert will be ligated into the vector as opposed to multiple insertions that can occur with blunt-ended cloning. In addition, with T-vector cloning there is no need to dephosphorylate the vector, and there is a low background of religated vector.

Using this procedure with optimized insert:vector ratios, 55–95% recombinants were obtained when *Pfu* and *Tli* DNA Polymerases were used to generate the insert DNA (Table 2). It is critical that the PCR fragments are purified using the Wizard® SV Gel and PCR Clean-Up System (Cat.# A9281) or by direct isolation from a gel by other means. In the absence of purification, the proofreading activity of the *Pfu*, *Pwo* and *Tli* DNA Polymerases will degrade the PCR fragments or remove the 3'-terminal deoxyadenosine added during tailing or the 3'-terminal deoxythymidine from the vector during the A-tailing reaction or ligation.

To optimize cloning efficiency, the amount of DNA in the A-tailing reaction and the ligation volumes must be adjusted depending on the molar yield of the purified PCR product. When molar concentrations are high due to small fragment size and/or good amplification, small volumes of the PCR fragment are needed for the A-tailing and ligation reactions. However, when molar concentration is low due to large fragment size and/or poor amplification, large volumes of the PCR fragment are needed for the A-tailing and ligation reactions. We have successfully used 1–7µl of the purified PCR fragment in the A-tailing reaction using *Taq* DNA Polymerase^(e) to optimize the insert:vector ratio. See Section IV.C for further discussion of optimizing the insert:vector ratio. Recombinants were identified by blue/white screening, and 70–100% of the recombinants were shown to have the correct size insert by PCR amplification of DNA. Few recombinants were observed in the control reactions in which the PCR fragment was not tailed. These control results confirm that the majority of the pGEM®-T Easy Vector used contains 3'-terminal deoxythymidine and that, during the A-tailing, *Taq* DNA Polymerase added a 3'-terminal deoxyadenosine to a significant proportion of the PCR fragment.

Table 2. Comparison of A-Tailing Procedures Used With Different DNA Polymerases.

Polymerase	% Recombinants ¹			
	1-Hour Ligation at 24°C (Standard)		16-Hour Ligation at 4°C (Alternative)	
	542bp	1.8kb	542bp	1.8kb
<i>Pfu</i> DNA Polymerase	65–84% ²	31–55% ³	81–95% ²	50–75% ³
<i>Tli</i> DNA Polymerase	68–77% ⁴	37–65% ⁵	85–93% ⁴	60–81% ⁵

PCR fragments generated by *Pfu* and *Tli* DNA Polymerase were A-tailed and ligated into pGEM®-T Easy Vector for 1 hour at 24°C or for 16 hours at 4°C. Two microliters of ligation mix was transformed into JM109 Competent Cells and plated on LB/amp/IPTG/X-gal plates.


¹% Recombinants = % white and/or pale blue colonies. PCR fragments were purified with the Wizard® PCR Preps DNA Purification System^(f) prior to A-tailing.

²Insert:vector ratios tested: 5:1, 3:1, 1:1. Volume of PCR amplification product used in A-tailing: 1–2µl.

³Insert:vector ratios tested: 3:1, 2:1, 1:1. Volume of PCR amplification product used in A-tailing: 3–7µl.

⁴Insert:vector ratios tested: 3:1, 2:1, 1:1. Volume of PCR amplification product used in A-tailing: 1–2µl.

⁵Insert:vector ratios tested: 2:1, 1:1. Volume of PCR amplification product used in A-tailing: 4–7µl.



2357MA02/BA

Figure 4. An A-tailing procedure for blunt-ended PCR fragments purified with Wizard® SV Gel and PCR Clean-Up System (Cat.# A9281) and used in T-vector cloning.

C. Optimizing Insert:Vector Molar Ratios

The pGEM®-T and pGEM®-T Easy Vector Systems have been optimized using a 1:1 molar ratio of the Control Insert DNA to the vectors. However, ratios of 8:1 to 1:8 have been used successfully. If initial experiments with your PCR product are suboptimal, ratio optimization may be necessary. Ratios from 3:1 to 1:3 provide good initial parameters. The concentration of PCR product should be estimated by comparison to DNA mass standards on a gel or by using a fluorescent assay (7). The pGEM®-T and pGEM®-T Easy Vectors are approximately 3kb and are supplied at 50ng/µl. To calculate the appropriate amount of PCR product (insert) to include in the ligation reaction, use the following equation.

$$\frac{\text{ng of vector} \times \text{kb size of insert}}{\text{kb size of vector}} \times \text{insert:vector molar ratio} = \text{ng of insert}$$

Example of insert:vector ratio calculation:

How much 0.5kb PCR product should be added to a ligation in which 50ng of 3.0kb vector will be used if a 3:1 insert:vector molar ratio is desired?

$$\frac{50\text{ng vector} \times 0.5\text{kb insert}}{3.0\text{kb vector}} \times \frac{3}{1} = 25\text{ng insert}$$

Note: Using the same parameters for a 1:1 insert:vector molar ratio, 8.3ng of a 0.5kb insert would be required.

D. Screening Transformants for Inserts

Successful cloning of an insert in the pGEM[®]-T and pGEM[®]-T Easy Vectors interrupts the coding sequence of β -galactosidase; recombinant clones can usually be identified by color screening on indicator plates. However, the characteristics of PCR products cloned into the pGEM[®]-T and pGEM[®]-T Easy Vectors can significantly affect the ratio of blue:white colonies obtained following transformation of competent cells. Clones that contain PCR products, in most cases, produce white colonies, but blue colonies can result from PCR fragments that are cloned in-frame with the *lacZ* gene. Such fragments are usually a multiple of 3 base pairs long (including the 3'-A overhangs), and do not contain in-frame stop codons. There have been reports of DNA fragments of up to 2kb that have been cloned in-frame and have produced blue colonies.

Even if your PCR product is not a multiple of 3 bases long, the amplification process can introduce mutations (e.g., deletions or point mutations) that may result in blue colonies when competent cells are transformed with the fragment inserted into the pGEM[®]-T or pGEM[®]-T Easy Vectors.

The Control Insert DNA supplied with the pGEM[®]-T and pGEM[®]-T Easy Systems is a 542bp fragment from pGEM[®]-*luc* Vector^(b,g) DNA (Cat.# E1541). This sequence has been mutated to contain multiple stop codons in all six reading frames, which ensures a low background of blue colonies for the control reaction. Results obtained with the Control Insert DNA may not be representative of those achieved with your PCR product.

E. Experimental Controls

Promega strongly recommends performing the controls detailed below. These are necessary to accurately assess the performance of the pGEM[®]-T and pGEM[®]-T Easy Vector Systems.

Positive Control

Set up a ligation reaction with the Control Insert DNA as described in the protocol (Section IV) and use it for transformations as described in Section V. This control will allow you to determine whether the ligation is proceeding efficiently. Typically, approximately 100 colonies should be observed, 10–40% of which are blue, when competent cells that have a transformation efficiency of 1×10^8 cfu/ μ g DNA are transformed. Greater than 60% of the colonies should be white. The Control Insert DNA is specifically designed to produce white colonies; however, other insert DNA may not yield white colonies (see Section VI.D). Background blue colonies from the positive control ligation reaction arise from non-T-tailed or undigested pGEM[®]-T or pGEM[®]-T Easy Vector. These blue colonies are a useful internal transformation control; if no colonies are obtained, the transformation has failed. If blue colonies are obtained, but no whites, the result suggests that the ligation reaction failed. If <50% white colonies are seen in this positive control reaction, then the ligation conditions were probably suboptimal.

The concentration of the Control Insert DNA is such that 2 μ l (4ng/ μ l) can be used in a 10 μ l ligation reaction to achieve a 1:1 molar ratio with 50ng of the pGEM[®]-T or pGEM[®]-T Easy Vectors.

Background Control

Set up a ligation reaction with 50ng of pGEM[®]-T or pGEM[®]-T Easy Vector and no insert as described in the protocol (Section IV.A) and use it for transformations as described in Section V. This ligation will allow determination of the number of background blue colonies resulting from non-T-tailed or undigested pGEM[®]-T or pGEM[®]-T Easy Vector alone. If the recommendations in Section V are followed closely, 10–30 blue colonies will typically be observed if the transformation efficiency of the competent cells is 1×10^8 cfu/ μ g DNA. (Under these conditions, cells that have an efficiency of 1×10^7 cfu/ μ g DNA would yield 1–3 blue colonies and cells with a transformation efficiency of 1×10^9 cfu/ μ g DNA would yield 100–300 blue colonies). Compare the number of blue colonies obtained with this background control to the number of blue colonies obtained in the standard reaction using the PCR product. If ligation of the PCR product yields dramatically more blue colonies than the background control reaction, then recombinants are probably among these blue colonies (see Section VI.D).

Transformation Control

Check the transformation efficiency of the competent cells by transforming them with an uncut plasmid (not pGEM[®]-T or pGEM[®]-T Easy since these vectors are linearized) and calculating cfu/ μ g DNA. If the transformation efficiency is lower than 1×10^8 cfu/ μ g DNA, prepare fresh cells. (Competent cells are available from Promega. See Section XII.B.) If you are not using JM109 High Efficiency Competent Cells (provided with pGEM[®]-T and pGEM[®]-T Easy Vector Systems II; Cat.# A3610 and A1380, respectively), be sure the cells are compatible with blue/white screening and standard ampicillin selection and have a transformation efficiency of at least 1×10^8 cfu/ μ g DNA.

Example of Transformation Efficiency Calculation:

After 100 μ l competent cells are transformed with 0.1ng uncut plasmid DNA, the transformation reaction is added to 900 μ l of SOC medium (0.1ng DNA/ml). From that volume, a 1:10 dilution with SOC medium (0.01ng DNA/ml) is made and 100 μ l plated on two plates (0.001ng DNA/100 μ l). If 200 colonies are obtained (average of two plates), what is the transformation efficiency?

$$\frac{200\text{cfu}}{0.001\text{ng}} = 2 \times 10^5\text{cfu/ng} = 2 \times 10^8\text{cfu}/\mu\text{g DNA}$$

VII. Isolation of Recombinant Plasmid DNA

A standard plasmid miniprep procedure, which takes 30–60 minutes to perform, is described in Promega's *Protocols and Applications Guide* (8). The miniprep process can be both laborious and time-consuming, particularly when large numbers of minipreps are required. A convenient and reliable method is the Wizard[®] Plus SV Minipreps DNA Purification System^(h,i) (Cat.# A1330).

VIII. Generation of Single-Stranded DNA from the pGEM[®]-T and pGEM[®]-T Easy Vectors

For induction of ssDNA production, bacterial cells containing either the pGEM[®]-T or pGEM[®]-T Easy Vector are infected with an appropriate helper phage (e.g., R408 Helper Phage, Cat.# P2291). The plasmid then enters the f1 replication mode, and the resulting ssDNA is exported as an encapsulated virus-like particle. The ssDNA is purified from the supernatant by simple precipitation and extraction procedures, which are described in detail in the *Protocols and Applications Guide* (Cat.# P1610). For further information, please contact your local Promega Branch Office or Distributor. In the U.S., contact Technical Services at 1-800-356-9526.

IX. Troubleshooting

Symptom	Possible Cause	Comments
No colonies	A problem has occurred with the transformation reaction or the cells have lost competence	Background undigested vector and religated non-T-tailed vector should yield 10–30 blue colonies independent of the presence of insert DNA. Check the background control (Section VI.E). Use high efficiency competent cells ($\geq 1 \times 10^8$ cfu/ μ g DNA). Test the efficiency by transforming the cells with an uncut plasmid that allows for antibiotic selection, such as the pGEM [®] -5Zf(+) Vector. If the guidelines in Section V.A are followed, cells at 1×10^8 cfu/ μ g DNA typically yield 100 colonies. Therefore, you would not see any colonies from cells that are $< 1 \times 10^7$ cfu/ μ g DNA (Section VI.E).
Less than 10% white colonies with Control Insert DNA	Improper dilution of the 2X Rapid Ligation Buffer	The T4 DNA ligase buffer is provided at a concentration of 2X. Use 5 μ l in a 10 μ l reaction.
	Ligation reaction has failed	Ligase buffer may have low activity. The 2X Rapid Ligation Buffer contains ATP, which degrades during temperature fluctuations. Avoid multiple freeze-thaw cycles by making single-use aliquots of the buffer. Use a fresh vial of buffer. To test the activity of the ligase and buffer set up a ligation with ~20ng of DNA markers (e.g., Lambda DNA/ <i>Hind</i> III Markers, Cat.# G1711). Compare ligated and nonligated DNA on a gel and check that the fragments have been religated into high molecular weight material.

IX. Troubleshooting (continued)

Symptom	Possible Cause	Comments
Less than 10% white colonies with Control Insert DNA (continued)	T-overhangs have been removed allowing blunt-ended ligation of vector and giving rise to more blue than white colonies	Avoid introduction of nucleases that may degrade the T-overhangs. Use only T4 DNA Ligase provided with the system, which has been tested for minimal exonuclease activity.
High colony number, but low percentage of white colonies with Control Insert DNA	Competent cells may have a high transformation efficiency ($\geq 1 \times 10^9$ cfu/ μ g) but there is a ligation problem	Approximately 1,000 colonies can be obtained in the positive control ligation using cells that are 10^9 cfu/ μ g DNA with 70–90% white colonies. If ligation is sub-optimal or fails, the total number of colonies will be high (up to 300 cells at 1×10^9 cfu/ μ g), but the amount of white colonies will be low or zero. See comments under “Ligation reaction has failed” (above).
Less than 60% white colonies with Control Insert DNA	Improper dilution of the 2X Rapid Ligation Buffer	The T4 DNA ligase buffer is provided at a concentration of 2X. Use 5 μ l in a 10 μ l reaction.
	T-overhangs have been removed allowing blunt-ended ligation of vector and giving rise to more blue than white colonies	Avoid introduction of nucleases that may degrade the T-overhangs. Use only T4 DNA Ligase provided with the system, which has been tested for minimal exonuclease activity.
	Ligation temperature is too high	Higher temperatures (>28°C) give rise to increased background and fewer recombinants.
Low number of or no white colonies with PCR product	Improper dilution of the 2X Rapid Ligation Buffer	The T4 DNA ligase buffer is provided at a concentration of 2X. Use 5 μ l in a 10 μ l reaction.
	Ligation incubation is not long enough	Optimal results are seen with an overnight ligation.
	Failed ligation due to an inhibitory component in the PCR product	Mix some of the PCR product with the positive control ligation to see if it is inhibiting the reaction. If an inhibitor is suspected, repurify the PCR fragment.
	PCR product is not ligating because there are no 3'-A overhangs	As summarized in Table 1, not all thermostable DNA polymerases create a 3'-A overhang (3,4). Blunt-ended fragments may be subsequently A-tailed by treatment with an appropriate polymerase and dATP (9,10).

For questions not addressed here, please contact your local Promega branch office or distributor (contact information available at: www.promega.com).

E-mail: techserv@promega.com

IX. Troubleshooting (continued)

Symptom	Possible Cause	Comments
Low number of or no white colonies with PCR product (continued)	PCR product cannot be ligated due to pyrimidine dimers formed from UV overexposure	This is a common problem with gel-purified DNA. There is no way to fix this; the DNA must be remade. Exposure to shortwave UV should be limited as much as possible. Use a glass plate between the gel and UV source to decrease UV overexposure. If possible, only visualize the PCR product using a longwave UV source.
	The PCR fragment is inserted, but it is not disrupting the <i>lacZ</i> gene	If there are a higher number of blue colonies resulting from the PCR fragment ligation than with the background control, some of these blue colonies may contain insert. Screen blue and pale blue colonies (see Section VI.C).
	Insert:vector ratio is not optimal	Check the integrity and quantity of your PCR fragment by gel analysis. Optimize the insert:vector ratio (see Section VI.C).
	There may be primer-dimers present in PCR fragment preparation	Primer-dimers will ligate into the pGEM [®] -T or pGEM [®] -T Easy Vector, but may not be seen after restriction digestion and gel analysis because of their small size. The vector will appear to contain no insert. More blue colonies may be seen with the ligation than on the background control plates. The PCR fragment should be gel-purified.
	Multiple PCR products are generated and cloned into the pGEM [®] -T or pGEM [®] -T Easy Vector.	Gel-purify the PCR fragment of interest.
	DNA has rearranged	Check a number of clones to see whether the rearrangement is random. If so, the clone of interest should be present and can be identified by screening several clones. If the same rearrangement is found in all of the clones, use a repair-deficient bacterial strain to protect the insert (e.g., SURE [®] cells), which may reduce recombination events.

IX. Troubleshooting (continued)

Symptom	Possible Cause	Comments
PCR product ligation reaction produces white colonies only (no blue colonies are present)	Ampicillin is inactive, allowing ampicillin-sensitive cells to grow	Check that ampicillin plates are made properly and used within 1 month. Test ampicillin activity by streaking plates, with and without ampicillin, using an ampicillin-sensitive clone.
	The bacterial strain (e.g., JM109) has lost its F' episome	Check the background control. If these colonies are not blue, the cells may have lost the F' episome (assuming <i>lacI^qZΔM15</i> is located on the F' in the transformed strain and appropriate plates were used). Be sure that the cells are prepared properly for use with this system (see Section V).
	Plates are incompatible with blue/white screening	Check the background control. If these colonies are not blue, check that the plates have ampicillin/IPTG/X-Gal and are fresh. If there is any question about the quality of the plates, repeat plating with fresh plates.
Not enough clones contain the PCR product of interest	Insufficient A-tailing of the PCR fragment	After purification of the PCR fragment, set up an A-tailing reaction (9,10). Clean up the sample and proceed with the protocol.
	Insert:vector ratio is not optimal	Check the integrity and quality of your PCR fragment by gel analysis. Optimize the insert:vector ratio (see Section VI.C).
	Multiple PCR products are generated and cloned into the pGEM [®] -T or pGEM [®] -T Easy Vector	Gel purify the PCR fragment of interest.

For questions not addressed here, please contact your local Promega branch office or distributor (contact information available at: www.promega.com).

E-mail: techserv@promega.com

X. References

1. Mezei, L.M. and Storts, D.R. (1994) Purification of PCR products. In: *PCR Technology: Current Innovations*, Griffin, H.G. and Griffin, A.M., eds., CRC Press, Boca Raton, FL, 21.
2. Robles, J. and Doers, M. (1994) pGEM[®]-T Vector Systems troubleshooting guide. *Promega Notes* **45**, 19–20.
3. Clark, J.M. (1988) Novel non-templated nucleotide addition reactions catalyzed by procaryotic and eucaryotic DNA polymerases. *Nucl. Acids Res.* **16**, 9677–86.
4. Newton, C.R. and Graham, A. (1994) In: *PCR*, BIOS Scientific Publishers, Ltd., Oxford, UK, 13.
5. Messing, J. *et al.* (1981) A system for shotgun DNA sequencing. *Nucl. Acids Res.* **9**, 309–21.
6. Knoche, K. and Kephart, D. (1999) Cloning blunt-end *Pfu* DNA Polymerase-generated PCR fragments into pGEM[®]-T Vector Systems. *Promega Notes* **71**, 10–13.
7. Haff, L. and Mezei, L. (1989) *Amplifications* **1**, 8.
8. *Protocols and Applications Guide*, Third Edition (1996) Promega Corporation.
9. Kobs, G. (1995) pGEM[®]-T Vector: cloning of modified blunt-ended DNA fragments. *Promega Notes* **55**, 28–29.
10. Zhou, M.-Y., Clark, S.E. and Gomez-Sanchez, C.E. (1995) Universal cloning method by TA strategy. *BioTechniques* **19**, 34–35.

XI. Appendix A: Vector Sequences and Restriction Sites

A. pGEM[®]-T Vector Sequence

The sequence supplied below is that of the circular pGEM[®]-5Zf(+) Vector from which the pGEM[®]-T Vector is derived. The pGEM[®]-T Vector has been linearized with *Eco*R V at base 51 of this sequence (indicated by an asterisk) and a T added to both 3'-ends. The added T is not included in this sequence. The sequence shown corresponds to RNA synthesized by T7 RNA Polymerase and is complementary to RNA synthesized by SP6 RNA Polymerase. The strand shown is complementary to the ssDNA produced by this vector. Vector sequences are also available at: www.promega.com/vectors/.

```

1   GGGCGAATTG  GGCCCGACGT  CGCATGCTCC  CGGCCGCCAT  GGCCGCGGGA
51  T*ATCACTAGT GCGGCCGCCT  GCAGGTCGAC  CATATGGGAG  AGCTCCCAAC
101 GCGTTGGATG  CATAGCTTGA  GTATTCTATA  GTGTCACCTA  AATAGCTTGG
151 CGTAATCATG  GTCATAGCTG  TTTCTGTGT  GAAATTGTTA  TCCGCTCACA
201 ATTCCACACA  ACATACGAGC  CGGAAGCATA  AAGTGTAAG  CCTGGGGTGC
251 CTAATGAGTG  AGCTAACTCA  CATTAATTGC  GTTGCCTCA  CTGCCCGCTT
301 TCCAGTCGGG  AAACCTGTCT  TGCCAGCTGC  ATTAATGAAT  CGGCCAACGC
351 GCGGGGAGAG  GCGGTTTTCG  TATTGGGCGC  TCTTCCGCTT  CCTCGCTCAC
401 TGACTCGCTG  CGCTCGGTCG  TTCGGCTGCG  GCGAGCGGTA  TCAGCTCACT
451 CAAAGGCGGT  AATACGGTTA  TCCACAGAAT  CAGGGGATAA  CGCAGGAAAG
501 AACATGTGAG  CAAAAGGCCA  GCAAAGGCC  AGGAACCGTA  AAAAGGCCGC
551 GTTGCTGGCG  TTTTTCATA  GGCTCCGCC  CCCTGACGAG  CATCACAAA
601 ATCGACGCTC  AAGTCAGAGG  TGGCGAAACC  CGACAGGACT  ATAAAGATAC

```

651	CAGGCGTTTC	CCCCTGGAAG	CTCCCTCGTG	CGCTCTCCTG	TTCCGACCCT
701	GCCGCTTACC	GGATACCTGT	CCGCCTTTCT	CCCTTCGGGA	AGCGTGCGC
751	TTTCTCATAG	CTCACGCTGT	AGGTATCTCA	GTTTCGGTGTA	GGTCGTTTCG
801	TCCAAGCTGG	GCTGTGTGCA	CGAACCCCCC	GTTTCAGCCCG	ACCGCTGCGC
851	CTTATCCGGT	AACTATCGTC	TTGAGTCCAA	CCCGGTAAGA	CACGACTTAT
901	CGCCACTGGC	AGCAGCCACT	GGTAACAGGA	TTAGCAGAGC	GAGGTATGTA
951	GGCGGTGCTA	CAGAGTTCTT	GAAGTGGTGG	CCTAACTACG	GCTACACTAG
1001	AAGAACAGTA	TTTGGTATCT	GCGCTCTGCT	GAAGCCAGTT	ACCTTCGGAA
1051	AAAGAGTTGG	TAGCTCTTGA	TCCGGCAAAC	AAACCACCGC	TGGTAGCGGT
1101	GGTTTTTTTG	TTTGCAAGCA	GCAGATTACG	CGCAGAAAAA	AAGGATCTCA
1151	AGAAGATCCT	TTGATCTTTT	CTACGGGGTC	TGACGCTCAG	TGGAACGAAA
1201	ACTCACGTTA	AGGGATTTTG	GTCATGAGAT	TATCAAAAAG	GATCTTCACC
1251	TAGATCCTTT	TAAATTAATA	ATGAAGTTTT	AAATCAATCT	AAAGTATATA
1301	TGAGTAAACT	TGGTCTGACA	GTTACCAATG	CTTAATCAGT	GAGGCACCTA
1351	TCTCAGCGAT	CTGTCTATTT	CGTTCATCCA	TAGTTGCCTG	ACTCCCCGTC
1401	GTGTAGATAA	CTACGATACG	GGAGGGCTTA	CCATCTGGCC	CCAGTGCTGC
1451	AATGATACCG	CGAGACCCAC	GCTCACCGGC	TCCAGATTTA	TCAGCAATAA
1501	ACCAGCCAGC	CGGAAGGGCC	GAGCGCAGAA	GTGGTCCTGC	AACTTTATCC
1551	GCCTCCATCC	AGTCTATTAA	TTGTTGCCGG	GAAGCTAGAG	TAAGTAGTTC
1601	GCCAGTTAAT	AGTTTGCGCA	ACGTTGTTGC	CATTGCTACA	GGCATCGTGG
1651	TGTCACGCTC	GTCGTTTGGT	ATGGCTTCAT	TCAGTCCCGG	TTCCCAACGA
1701	TCAAGGCGAG	TTACATGATC	CCCCATGTTG	TGCAAAAAAG	CGGTTAGCTC
1751	CTTCGGTCCT	CCGATCGTTG	TCAGAAGTAA	GTTGGCCGCA	GTGTTATCAC
1801	TCATGGTTAT	GGCAGCACTG	CATAATTCTC	TTACTGTCAT	GCCATCCGTA
1851	AGATGCTTTT	CTGTGACTGG	TGAGTACTCA	ACCAAGTCAT	TCTGAGAATA
1901	GTGTATGCGG	CGACCGAGTT	GCTCTTGCCC	GGCGTCAATA	CGGGATAATA
1951	CCGCGCCACA	TAGCAGAACT	TTAAAAGTGC	TCATCATTGG	AAAACGTTCT
2001	TCGGGGCGAA	AACTCTCAAG	GATCTTACCG	CTGTTGAGAT	CCAGTTCGAT
2051	GTAACCCACT	CGTGCACCCA	ACTGATCTTC	AGCATCTTTT	ACTTTCACCA
2101	GCGTTTCTGG	GTGAGCAAAA	ACAGGAAGGC	AAAATGCCGC	AAAAAAGGGA
2151	ATAAGGGCGA	CACGGAAATG	TTGAATACTC	ATACTCTTCC	TTTTTCAATA
2201	TTATTGAAGC	ATTTATCAGG	GTTATTGTCT	CATGAGCGGA	TACATATTTG
2251	AATGTATTTA	GAAAAATAAA	CAAATAGGGG	TTCCGCGCAC	ATTTCCCCGA
2301	AAAGTGCCAC	CTGATGCGGT	GTGAAATACC	GCACAGATGC	GTAAGGAGAA
2351	AATACCGCAT	CAGGAAATTG	TAAGCGTTAA	TATTTTGTTA	AAATTCGCGT
2401	TAAATTTTTG	TTAAATCAGC	TCATTTTTTA	ACCAATAGGC	CGAAATCGGC
2451	AAAATCCCTT	ATAAATCAAA	AGAATAGACC	GAGATAGGGT	TGAGTGTGTT
2501	TCCAGTTTGG	AACAAGAGTC	CACTATTAAT	GAACGTGGAC	TCCAACGTCA
2551	AAGGGCGAAA	AACCGTCTAT	CAGGGCGATG	GCCCACTACG	TGAACCATCA

```

2601 CCCTAATCAA GTTTTTTGGG GTCGAGGTGC CGTAAAGCAC TAAATCGGAA
2651 CCCTAAAGGG AGCCCCGAT TTAGAGCTTG ACGGGGAAAG CCGGCGAACG
2701 TGGCGAGAAA GGAAGGGAAG AAAGCGAAAG GAGCGGGCGC TAGGGCGCTG
2751 GCAAGTGTAG CGGTCACGCT GCGCGTAACC ACCACACCCG CCGCGCTTAA
2801 TGCGCCGCTA CAGGGCGCGT CCATTGCGCA TTCAGGCTGC GCAACTGTTG
2851 GGAAGGGCGA TCGGTGCGGG CCTCTTCGCT ATTACGCCAG CTGGCGAAAG
2901 GGGGATGTGC TGCAAGGCGA TTAAGTTGGG TAACGCCAGG GTTTTCCCAG
2951 TCACGACGTT GTAAAACGAC GGCCAGTGAA TTGTAATACG ACTCACTATA

```

B. pGEM®-T Vector Restriction Sites

The following restriction enzyme tables are based on those of the circular pGEM®-5Zf(+) Vector from which the pGEM®-T Vector is derived. The pGEM®-T Vector has been linearized at base 51 with *EcoR V* and a T added to both 3'-ends. This site will not be recovered upon ligation of the vector and insert. The tables were constructed using DNASTAR® sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3'-end of the cut DNA (the base to the left of the cut site). Please contact your local Promega Branch Office or Distributor if you identify a discrepancy. In the U.S., contact Technical Services at 1-800-356-9526.

Table 3. Restriction Enzymes That Cut the pGEM®-T Vector Between 1 and 5 Times.

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Aat II	1	20	Bst Z I	2	31, 62
Acc I	1	76	<i>Cfr10 I</i>	2	1475, 2690
Acy I	2	17, 1932	Dde I	4	777, 1186, 1352, 1892
<i>Afl III</i>	2	99, 502	Dra I	3	1261, 1280, 1972
Alw26 I	2	1456, 2232	<i>Dra III</i>	1	2589
Alw44 I	2	816, 2062	<i>Drd I</i>	2	610, 2544
<i>AlwN I</i>	1	918	<i>Dsa I</i>	2	37, 43
Apa I	1	14	<i>Eag I</i>	2	31, 62
<i>AspH I</i>	4	94, 820, 1981, 2066	<i>Ear I</i>	3	386, 2190, 2878
Ava II	2	1533, 1755	EclHK I	1	1395
Ban I	3	246, 1343, 2626	Eco52 I	2	31, 62
Ban II	3	14, 94, 2664	EcoCR I	1	92
Bbu I	1	26	EcoR V	1	51 (see above)
Bgl I	3	39, 1515, 2833	Fok I	5	119, 1361, 1542, 1829, 2919
<i>Bsa I</i>	1	1456	<i>Fsp I</i>	2	1617, 2840
<i>BsaA I</i>	1	2589	Hae II	4	380, 750, 2740, 2748
<i>BsaH I</i>	2	17, 1932	<i>Hga I</i>	4	613, 1191, 1921, 2806
<i>BsaJ I</i>	5	37, 43, 241, 662, 2936	Hinc II	1	77
<i>Bsp120 I</i>	1	10	<i>Hind II</i>	1	77
<i>BspH I</i>	2	1222, 2230	Hsp92 I	2	17, 1932
<i>BspM I</i>	1	62	<i>Mae I</i>	5	56, 997, 1250, 1585, 2740
<i>BssS I</i>	2	675, 2059	Mlu I	1	99
BstO I	5	242, 530, 651, 664, 2937			
BstX I	1	103			

Note: The enzymes listed in boldface type are available from Promega.

Table 3. Restriction Enzymes That Cut the pGEM®-T Vector Between 1 and 5 Times (continued).

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Nae I	1	2692	Sac II	1	46
Nci I	4	30, 882, 1578, 1929	Sal I	1	75
Nco I	1	37	Sca I	1	1875
Nde I	1	82	Sfi I	1	39
NgoM IV	1	2690	Sin I	2	1533, 1755
Not I	1	62	Spe I	1	55
Nsi I	1	112	Sph I	1	26
<i>Nsp I</i>	2	26, 506	<i>Sse8387 I</i>	1	73
<i>Ppu10 I</i>	1	108	Ssp I	2	2199, 2381
Pst I	1	73	Sty I	1	37
Pvu I	2	1765, 2861	Taq I	4	76, 602, 2046, 2622
Pvu II	2	326, 2890	<i>Tfi I</i>	2	337, 477
Rsa I	1	1875	Vsp I	3	273, 332, 1567
Sac I	1	94	Xmn I	1	1994

Note: The enzymes listed in boldface type are available from Promega.

Table 4. Restriction Enzymes That Do Not Cut the pGEM®-T Vector.

AccB7 I	<i>Bbs I</i>	Bst98 I	<i>Ehe I</i>	<i>PflM I</i>	SnaB I
Acc III	Bcl I	BstE II	<i>Fse I</i>	<i>PinA I</i>	<i>SpI I</i>
Acc65 I	Bgl II	Bsu36 I	Hind III	<i>Pme I</i>	<i>Srf I</i>
<i>Afl II</i>	<i>Blp I</i>	Cla I	Hpa I	<i>Pml I</i>	Stu I
Age I	<i>Bpu1102 I</i>	Csp I	I-Ppo I	<i>PpuM I</i>	<i>Swa I</i>
<i>Asc I</i>	<i>BsaB I</i>	Csp45 I	<i>Kas I</i>	<i>PshA I</i>	Tth111 I
Ava I	BsaM I	<i>Dra II</i>	Kpn I⁽ⁱ⁾	<i>Psp5 II</i>	Xba I
<i>Avr II</i>	<i>Bsm I</i>	Eco47 III	Nar I	<i>PspA I</i>	<i>Xcm I</i>
Bal I	BsrBR I	<i>Eco72 I</i>	Nhe I	<i>Rsr II</i>	Xho I
BamH I	<i>BsrG I</i>	<i>Eco81 I</i>	Nru I	Sgf I^(k)	Xma I
<i>Bbe I</i>	BssH II	<i>EcoN I</i>	<i>Pac I</i>	<i>SgrA I</i>	
<i>BbrP I</i>	<i>Bst1107 I</i>	EcoR I	<i>PaeR7 I</i>	Sma I	

Note: The enzymes listed in boldface type are available from Promega.

Table 5. Restriction Enzymes That Cut the pGEM®-T Vector 6 or More Times.

<i>Aci I</i>	Bst71 I	Hae III	<i>Mae III</i>	Nde II	<i>SfaN I</i>
Alu I	<i>BstU I</i>	Hha I	Mbo I	<i>Nla III</i>	Tru9 I
<i>Bbv I</i>	Cfo I	Hinf I	Mbo II	<i>Nla IV</i>	Xho II
BsaO I	Dpn I	Hpa II	<i>Mnl I</i>	<i>Ple I</i>	
Bsp1286 I	<i>Dpn II</i>	<i>Hph I</i>	<i>Mse I</i>	Sau3A I	
<i>Bsr I</i>	<i>Eae I</i>	Hsp92 II	Msp I	Sau96 I	
BsrS I	<i>Fnu4H I</i>	<i>Mae II</i>	MspA1 I	<i>ScrF I</i>	

Note: The enzymes listed in boldface type are available from Promega.

C. pGEM[®]-T Easy Vector Sequence

The pGEM[®]-T Easy Vector has been linearized with *EcoR* V at base 60 of this sequence (indicated by an asterisk) and a T added to both 3'-ends. The added T is not included in this sequence. The sequence shown corresponds to RNA synthesized by T7 RNA Polymerase and is complementary to RNA synthesized by SP6 RNA Polymerase. The strand shown is complementary to the ssDNA produced by this vector. Vector sequences are also available at: www.promega.com/vectors/.

```

1  GGGCGAATTG  GGCCCGACGT  CGCATGCTCC  CGGCCGCCAT  GGCGGCCGCG
51  GGAATTCGAT* ATCACTAGTG  AATTCGCGGC  CGCCTGCAGG  TCGACCATAT
101 GGGAGAGCTC  CCAACGCGTT  GGATGCATAG  CTTGAGTATT  CTATAGTGTC
151 ACCTAAATAG  CTTGGCGTAA  TCATGGTCAT  AGCTGTTTCC  TGTGTGAAAT
201 TGTTATCCGC  TCACAATTCC  ACACAACATA  CGAGCCGGAA  GCATAAAGTG
251 TAAAGCCTGG  GGTGCCTAAT  GAGTGAGCTA  ACTCACATTA  ATTGCGTTGC
301 GCTCACTGCC  CGCTTTCAG  TCGGGAAACC  TGTCGTGCCA  GCTGCATTAA
351 TGAATCGGCC  AACGCGCGGG  GAGAGGCGGT  TTGCGTATTG  GGCGCTCTTC
401 CGCTTCCTCG  CTCCTGACT  CGCTGCGCTC  GGTGTTCCGG  CTGCGGCGAG
451 CGGTATCAGC  TCACTCAAAG  GCGGTAATAC  GGTATCCAC  AGAATCAGGG
501 GATAACGCAG  GAAAGAACAT  GTGAGCAAAA  GGCCAGCAA  AGGCCAGGAA
551 CCGTAAAAAG  GCCGCGTTGC  TGGCGTTTTT  CCATAGGCTC  CGCCCCCTG
601 ACGAGCATCA  CAAAAATCGA  CGCTCAAGTC  AGAGGTGGCG  AAACCCGACA
651 GGACTATAAA  GATACCAGGC  GTTTCCCCCT  GGAAGCTCCC  TCGTGCGCTC
701 TCCTGTTCCG  ACCCTGCCGC  TTACCGGATA  CCTGTCCGCC  TTTCTCCCTT
751 CGGGAAGCGT  GGCCTTTCT  CATAGCTCAC  GCTGTAGGTA  TCTCAGTTCG
801 GTGTAGGTCG  TTCGCTCAA  GCTGGGCTGT  GTGCACGAAC  CCCCCGTTCA
851 GCCCGACCGC  TCGCCTTAT  CCGGTAACTA  TCGTCTTGAG  TCCAACCCGG
901 TAAGACACGA  CTTATCGCCA  CTGGCAGCAG  CCACTGGTAA  CAGGATTAGC
951 AGAGCGAGGT  ATGTAGGCGG  TGCTACAGAG  TTCTTGAAGT  GGTGGCCTAA
1001 CTACGGCTAC  ACTAGAAGAA  CAGTATTTGG  TATCTGCGCT  CTGCTGAAGC
1051 CAGTTACCTT  CGGAAAAGA  GTTGGTAGCT  CTTGATCCGG  CAAACAAACC
1101 ACCGCTGGTA  GCGGTGGTTT  TTTTGTTCG  AAGCAGCAGA  TTACGCGCAG
1151 AAAAAAAGGA  TCTCAAGAAG  ATCCTTTGAT  CTTTTCTACG  GGGTCTGACG
1201 CTCAGTGGAA  CGAAACTCA  CGTTAAGGGA  TTTTGGTCAT  GAGATTATCA
1251 AAAAGGATCT  TCACCTAGAT  CCTTTTAAAT  TAAAAATGAA  GTTTTAAATC
1301 AATCTAAAGT  ATATATGAGT  AACTTTGGTC  TGACAGTTAC  CAATGCTTAA
1351 TCAGTGAGGC  ACCTATCTCA  GCGATCTGTC  TATTTTCGTT  ATCCATAGTT
1401 GCCTGACTCC  CCGTCGTGTA  GATAACTACG  ATACGGGAGG  GCTTACCATC
1451 TGGCCCCAGT  GCTGCAATGA  TACCGCGAGA  CCCACGCTCA  CCGGCTCCAG
1501 ATTTATCAGC  AATAAACAG  CCAGCCGGAA  GGGCCGAGCG  CAGAAGTGGT
1551 CCTGCAACTT  TATCCGCCTC  CATCCAGTCT  ATTAATTGTT  GCCGGGAAGC
1601 TAGAGTAAGT  AGTTCGCCAG  TTAATAGTTT  GCGCAACGTT  GTTGCCATTG

```

```

1651 CTACAGGCAT CGTGGTGTCA CGCTCGTCGT TTGGTATGGC TTCATTCAGC
1701 TCCGGTCCC AACGATCAAG GCGAGTTACA TGATCCCCCA TGTTGTGCAA
1751 AAAAGCGGTT AGCTCCTTCG GTCCTCCGAT CGTTGTCAGA AGTAAGTTGG
1801 CCGCAGTGTT ATCACTCATG GTTATGGCAG CACTGCATAA TTCTCTTACT
1851 GTCATGCCAT CCGTAAGATG CTTTTCTGTG ACTGGTGAGT ACTCAACCAA
1901 GTCATTCTGA GAATAGTGTA TGCGGCGACC GAGTTGCTCT TGCCCCGCGT
1951 CAATACGGGA TAATACCGCG CCACATAGCA GAACTTTAAA AGTGCTCATC
2001 ATTGGAAAAC GTTCTTCGGG GCGAAAATC TCAAGGATCT TACCGCTGTT
2051 GAGATCCAGT TCGATGTAAC CCACTCGTGC ACCCAACTGA TCTTCAGCAT
2101 CTTTTACTTT CACCAGCGTT TCTGGGTGAG CAAAAACAGG AAGGCAAAAT
2151 GCCGCAAAA AGGAATAAG GCGACACGG AAATGTTGAA TACTCATACT
2201 CTTCTTTTT CAATATTATT GAAGCATTTA TCAGGGTTAT TGTCTCATGA
2251 GCGGATACAT ATTTGAATGT ATTTAGAAAA ATAAACAAAT AGGGGTCCG
2301 CGCACATTT CCCGAAAAGT GCCACCTGAT GCGGTGTGAA ATACCGCACA
2351 GATGCGTAAG GAGAAAATAC CGCATCAGGA AATTGTAAGC GTTAATATTT
2401 TGTTAAAATT CGCGTTAAAT TTTTGTAAA TCAGCTCATT TTTTAACCAA
2451 TAGGCCGAAA TCGGCAAAAT CCCTTATAAA TCAAAAGAAT AGACCGAGAT
2501 AGGGTTGAGT GTTGTTCCAG TTTGGAACAA GAGTCCACTA TTAAAGAACG
2551 TGGACTCAA CGTCAAAGGG CGAAAAACCG TCTATCAGGG CGATGGCCCA
2601 CTACGTGAAC CATCACCTA ATCAAGTTTT TTGGGGTCGA GGTGCCGTAA
2651 AGCACTAAAT CGGAACCCTA AAGGGAGCCC CCGATTTAGA GCTTGACGGG
2701 GAAAGCCGGC GAACGTGGCG AGAAAGGAAG GGAAGAAAGC GAAAGGAGCG
2751 GGCGCTAGGG CGCTGGCAAG TGTAGCGGTC ACGCTGCGCG TAACCACCAC
2801 ACCCGCCGCG CTTAATGCGC CGCTACAGGG CGCGTCCATT CGCCATTCAG
2851 GCTGCGCAAC TGTTGGGAAG GCGGATCGGT GCGGGCCTCT TCGCTATTAC
2901 GCCAGCTGGC GAAAGGGGGA TGTGCTGCAA GCGGATTAAG TTGGGTAAACG
2951 CCAGGGTTTT CCCAGTCACG ACGTTGTAAA ACGACGGCCA GTGAATTGTA
3001 ATACGACTCA CTATA

```

D. pGEM[®]-T Easy Vector Restriction Sites

The pGEM[®]-T Easy Vector has been linearized at base 60 with *EcoR V* and a T added to both 3'-ends. This site will not be recovered upon ligation of the vector and insert. The tables were constructed using DNASTAR[®] sequence analysis software. Please note that we have not verified this information by restriction digestion with each enzyme listed. The location given specifies the 3'-end of the cut DNA (the base to the left of the cut site). Please contact your local Promega Branch Office or Distributor if you identify a discrepancy. In the U.S., contact Technical Services at 1-800-356-9526.

Table 6. Restriction Enzymes That Cut the pGEM[®]-T Easy Vector Between 1 and 5 Times

Enzyme	# of Sites	Location	Enzyme	# of Sites	Location
Aat II	1	20	Fok I	5	134, 1376, 1557, 1844, 2931
Acc I	1	91	<i>Fsp I</i>	2	1632, 2855
Acy I	2	17, 1947	Hae II	4	395, 765, 2755, 2763
<i>Afl III</i>	2	114, 517	<i>Hga I</i>	4	628, 1206, 1936, 2821
Alw26 I	2	1471, 2247	Hinc II	1	92
Alw44 I	2	831, 2077	<i>Hind II</i>	1	92
<i>AlwNI</i>	1	933	Hsp92 I	2	17, 1947
Apa I	1	14	<i>Mae I</i>	5	65, 1012, 1265, 1600, 2755
<i>AspH I</i>	4	109, 835, 1996, 2081	Mlu I	1	114
Ava II	2	1548, 1770	Nae I	1	2707
Ban I	3	261, 1358, 2641	Nci I	4	30, 897, 1593, 1944
Ban II	3	14, 109, 2679	Nco I	1	37
Bbu I	1	26	Nde I	1	97
Bgl I	4	39, 42, 1530, 2848	NgoM IV	1	2705
<i>Bsa I</i>	1	1471	Not I	2	43, 77
<i>BsaA I</i>	1	2604	Nsi I	1	127
<i>BsaH I</i>	2	17, 1947	<i>Nsp I</i>	2	26, 521
<i>BsaJ I</i>	5	37, 46, 256, 677, 2951	<i>Ppu10 I</i>	1	123
<i>Bsp120 I</i>	1	10	Pst I	1	88
<i>BspH I</i>	2	1237, 2245	Pvu I	2	1780, 2876
<i>BspM I</i>	1	77	Pvu II	2	341, 2905
<i>BssS I</i>	2	690, 2074	Rsa I	1	1890
BstO I	5	257, 545, 666, 679, 2952	Sac I	1	109
BstX I	1	118	Sac II	1	49
BstZ I	3	31, 43, 77	Sal I	1	90
<i>Cfr10 I</i>	2	1490, 2705	Sca I	1	1890
Dde I	4	792, 1201, 1367, 1907	Sin I	2	1548, 1770
Dra I	3	1276, 1295, 1987	Spe I	1	64
<i>Dra III</i>	1	2604	Sph I	1	26
<i>Drd I</i>	2	625, 2559	<i>Sse8387 I</i>	1	88
<i>Dsa I</i>	2	37, 46	Ssp I	2	2214, 2396
<i>Eag I</i>	3	31, 43, 77	Sty I	1	37
<i>Ear I</i>	3	401, 2205, 2893	Taq I	5	56, 91, 617, 2061, 2637
EcI/HK I	1	1410	<i>Tfi I</i>	2	352, 492
Eco52 I	3	31, 43, 77	Vsp I	3	288, 347, 1582
EcoICR I	1	107	Xmn I	1	2009
EcoR I	2	52, 70			
EcoR V	1	60 (see above)			

Note: The enzymes listed in boldface type are available from Promega.

Table 7. Restriction Enzymes That Do Not Cut the pGEM®-T Easy Vector.

Acc B7 I	<i>Bbs I</i>	Bst98 I	<i>Fse I</i>	<i>PinA I</i>	SnaB I
Acc III	Bcl I	BstE II	Hind III	<i>Pme I</i>	<i>Spl I</i>
Acc65 I	Bgl II	Bsu36 I	Hpa I	<i>Pml I</i>	<i>Srf I</i>
<i>Afl II</i>	<i>Blp I</i>	Clal	I-Ppo I	<i>PpuM I</i>	Stu I
Age I	<i>Bpu 1102 I</i>	Csp I	<i>Kas I</i>	<i>PshA I</i>	<i>Swa I</i>
<i>Asc I</i>	<i>BsaB I</i>	Csp45 I	Kpn I	<i>Psp5 II</i>	Tth111 I
Ava I	BsaM I	<i>Dra II</i>	Nar I	<i>PspA I</i>	Xba I
<i>Avr II</i>	<i>Bsm I</i>	Eco47 III	Nhe I	<i>Rsr II</i>	<i>Xcm I</i>
Bal I	Bsr BR I	<i>Eco72 I</i>	Nru I	Sfi I	Xho I
BamH I	<i>BsrG I</i>	<i>Eco81 I</i>	<i>Pac I</i>	Sgf I^(k)	Xma I
<i>Bbe I</i>	BssH II	<i>EcoN I</i>	<i>PaeR7 I</i>	<i>SgrA I</i>	
<i>BbrP I</i>	<i>Bst1107 I</i>	<i>Ehe I</i>	<i>PflM I</i>	Sma I	

Note: The enzymes listed in boldface type are available from Promega.

Table 8. Restriction Enzymes That Cut the pGEM®-T Easy Vector 6 or More Times.

<i>Aci I</i>	Bst71 I	Hae III	<i>Mae III</i>	Nde II	<i>SfaN I</i>
Alu I	<i>BstU I</i>	Hha I	Mbo I	<i>Nla III</i>	Tru9 I
<i>Bbv I</i>	Cfo I	Hinf I	Mbo II	<i>Nla IV</i>	Xho II
BsaO I	Dpn I	Hpa II	<i>Mnl I</i>	<i>Ple I</i>	
Bsp1286 I	<i>Dpn II</i>	<i>Hph I</i>	<i>Mse I</i>	Sau3A I	
<i>Bsr I</i>	<i>Eae I</i>	Hsp92 II	Msp I	Sau96 I	
BsrS I	<i>Fnu4H I</i>	<i>Mae II</i>	MspA1 I	<i>ScrF I</i>	

Note: The enzymes listed in boldface type are available from Promega.

XII. Appendix B: Reference Information

A. Composition of Buffers and Solutions

IPTG stock solution (0.1M)

1.2g IPTG (Cat.# V3951)

Add water to 50ml final volume. Filter-sterilize and store at 4°C.

X-Gal (2ml)

100mg 5-bromo-4-chloro-3-indolyl-β-D-galactoside (Cat.# V3941)

Dissolve in 2ml N,N'-dimethylformamide. Cover with aluminum foil and store at -20°C.

LB medium (per liter)

10g Bacto®-tryptone
5g Bacto®-yeast extract
5g NaCl

Adjust pH to 7.0 with NaOH.

LB plates with ampicillin

Add 15g agar to 1 liter of LB medium. Autoclave. Allow the medium to cool to 50°C before adding ampicillin to a final concentration of 100µg/ml. Pour 30–35ml of medium into 85mm petri dishes. Let the agar harden. Store at 4°C for up to 1 month or at room temperature for up to 1 week.

LB plates with ampicillin/IPTG/X-Gal

Make the LB plates with ampicillin as above; then supplement with 0.5mM IPTG and 80µg/ml X-Gal and pour the plates. Alternatively, 100µl of 100mM IPTG and 20µl of 50mg/ml X-Gal may be spread over the surface of an LB-ampicillin plate and allowed to absorb for 30 minutes at 37°C prior to use.

A. Composition of Buffers and Solutions (continued)

SOC medium (100ml)

2.0g	Bacto®-tryptone
0.5g	Bacto®-yeast extract
1ml	1M NaCl
0.25ml	1M KCl
1ml	2M Mg ²⁺ stock, filter-sterilized
1ml	2M glucose, filter-sterilized

Add Bacto®-tryptone, Bacto®-yeast extract, NaCl and KCl to 97ml distilled water. Stir to dissolve. Autoclave and cool to room temperature. Add 2M Mg²⁺ stock and 2M glucose, each to a final concentration of 20mM. Bring to 100ml with sterile, distilled water. The final pH should be 7.0.

2M Mg²⁺ stock

20.33g	MgCl ₂ • 6H ₂ O
24.65g	MgSO ₄ • 7H ₂ O

Add distilled water to 100ml. Filter sterilize.

2X Rapid Ligation Buffer, T4 DNA Ligase (provided)

60mM	Tris-HCl (pH 7.8)
20mM	MgCl ₂
20mM	DTT
2mM	ATP
10%	polyethylene glycol (MW8000, ACS Grade)

Store in single-use aliquots at –20°C. Avoid multiple freeze-thaw cycles.

TYP broth (per liter)

16g	Bacto®-tryptone
16g	Bacto®-yeast extract
5g	NaCl
2.5g	K ₂ HPO ₄

B. Related Products

PCR Cloning Systems

Product	Size	Cat.#
pTARGET™ Mammalian Expression Vector System ^(a,l)	20 reactions	A1410
Direct mammalian expression from a T-Vector.		

Product	Size	Cat.#
PinPoint™ Xa-1 T-Vector System I ^(a,b,m)	20 reactions	V2610
PinPoint™ Xa-1 T-Vector System II with Competent Cells ^(a,b,m)	20 reactions	V2850
Direct bacterial expression from a T-vector.		

Thermostable Polymerases

Product	Size	Cat.#
PCR Master Mix ^(n,o)	100 reactions	M7501
	1,000 reactions	M7502

Premixed 2X solution containing *Taq* DNA Polymerase, dNTPs, reaction buffer and magnesium. Simply add template and primers. For Laboratory Use.

Product	Size	Cat.#
GoTaq™ DNA Polymerase*(e,o)	100u	M3001
	500u	M3005
<i>Taq</i> DNA Polymerase*(e)	100u	M1661
	500u	M1665
<i>Pfu</i> DNA Polymerase(e) (not available in the U.S.)	100u	M7741
	500u	M7745
<i>Tli</i> DNA Polymerase*(e)	50u	M7101

*For Laboratory Use.

RT-PCR Systems

Product	Size	Cat.#
Access RT-PCR System ⁽ⁿ⁾	20 reactions	A1260
	100 reactions	A1250
	500 reactions	A1280
AccessQuick™ RT-PCR System ^(n,p)	20 reactions	A1701
	100 reactions	A1702
	500 reactions	A1703
Improm-II™ Reverse Transcription System	100 reactions	A3800

For Laboratory Use.

PCR Purification Systems

Product	Size	Cat.#
Wizard® SV Gel and PCR Clean-Up System ^(d)	50 preps	A9281
	250 preps	A9282

For Laboratory Use.

Product	Size	Cat.#
Wizard® SV 96 PCR Clean-Up System	1 × 96 preps	A9340
	4 × 96 preps	A9341
	8 × 96 preps	A9342
Wizard® MagneSil™ PCR Clean-Up System ^(q)	4 × 96 preps	A1930
	8 × 96 preps	A1931
	100 × 96 preps	A1935

For Laboratory Use.

dNTPs

Product	Size	Cat.#
PCR Nucleotide Mix (10mM each)	200µl	C1141
	1,000µl	C1145
dATP, dCTP, dGTP, dTTP, each at 100mM	10µmol of each	U1330
dATP, dCTP, dGTP, dTTP, each at 100mM	40µmol of each	U1240
dATP, dCTP, dGTP, dTTP, each at 100mM	200µmol of each	U1410

For Laboratory Use.

Sequencing Primers

Product	Size	Cat.#
SP6 Promoter Primer	2µg	Q5011
T7 Promoter Primer	2µg	Q5021
pUC/M13 Primer, Forward (24mer)	2µg	Q5601
pUC/M13 Primer, Reverse (22mer)	2µg	Q5421

- (a) Licensed under one or more of U.S. Pat. No. 5,487,993 and European Pat. No. 0 550 693.
- (b) U.S. Pat. No. 4,766,072.
- (c) U.S. Pat. No. 5,552,302, Australian Pat. No. 646803 and other patents. U.S. Pat. Nos. 4,966,964, 5,019,556 and 5,266,687, Australian Pat. Nos. 61688, 641261 and other pending and issued patents, which claim vectors encoding a portion of human placental ribonuclease inhibitor, are exclusively licensed to Promega Corporation.
- (d) U.S. Pat. Nos. 5,658,548, 5,808,041, Australian Pat. No. 689815 and other patents pending.
- (e) Certain applications of this product are covered by patents issued and applicable in certain countries. Because purchase of this product does not include a license to perform any patented application, users of this product may be required to obtain a patent license depending upon the particular application and country in which the product is used.
- (f) Licensed under U.S. Pat. No. 5,075,430.
- (g) The method of recombinant expression of *Coleoptera* luciferase is covered by U.S. Pat. Nos. 5,583,024, 5,674,713 and 5,700,673.
- (h) Australian Pat. No. 730718 and other patents and patents pending.
- (i) U.S. Pat. No. 5,981,235, Australian Pat. No. 729932 and other patents pending.
- (j) Purchase of this product is accompanied by a limited license under U.S. Pat. Nos. 5,082,784 and 5,192,675 for the internal research use of the purchaser.
- (k) U.S. Pat. No. 5,391,487.
- (l) The CMV promoter and its use are covered under U.S. Pat. Nos. 5,168,062 and 5,385,839 owned by the University of Iowa Research Foundation, Iowa City, Iowa, and licensed FOR RESEARCH USE ONLY. Commercial users must obtain a license to these patents directly from the University of Iowa Research Foundation.
- (m) For research purposes only. Not for diagnostic or therapeutic use. For nonresearch uses of the portion of the vector encoding the biotinylation sequence, please contact Promega Corporation for licensing information. For bulk purchases of the SoftLink™ Resin, contact TosoHaas, 156 Keystone Drive, Montgomeryville, PA 18936, 1-800-366-4875 or 215-283-5000.
- (n) The PCR process is covered by patents issued and applicable in certain countries. Promega does not encourage or support the unauthorized or unlicensed use of the PCR process. Use of this product is recommended for persons that either have a license to perform PCR or are not required to obtain a license.
- (o) U.S. Pat. No. 6,242,235 and other patents pending.
- (p) U.S. Pat. Nos. 4,966,964, 5,019,556 and 5,266,687, which claim vectors encoding a portion of human placental ribonuclease inhibitor, are exclusively licensed to Promega Corporation.
- (q) U.S. Pat. Nos. 6,027,945, 6,368,800, Australian Pat. No. 732756, and other patents and patents pending.

© 1996–2003 Promega Corporation. All Rights Reserved.

Erase-a-Base, pGEM, Riboprobe and Wizard are trademarks of Promega Corporation and are registered with the U.S. Patent and Trademark Office. AccessQuick, ImProm-II, MagneSil, PinPoint and pTARGET are trademarks of Promega Corporation.

AmpliQ is a registered trademark of Roche Molecular Systems, Inc. Bacto is a registered trademark of Difco Laboratories. DNASTAR is a registered trademark of DNASTAR, Inc. Luer-Lok is a registered trademark of Becton Dickinson and Company. SURE is a registered trademark of Stratagene Cloning Systems. Vent is a registered trademark of New England Biolabs, Inc.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.