



Platinum® Pfx DNA Polymerase

Cat. No. 11708-013	Size: 100 reactions
11708-021	250 reactions
11708-039	500 reactions

Conc: 2.5 U/μl	Store at -20°C
-----------------------	-----------------------

Note: Read and follow reaction conditions carefully to ensure optimal performance.

Description

Platinum® Pfx DNA Polymerase is a proprietary enzyme preparation containing recombinant DNA polymerase from *Thermococcus* sp. strain KOD (1,2). Platinum® Pfx DNA Polymerase possesses proofreading 3' to 5' exonuclease activity and provides higher fidelity than Pfu DNA polymerase (3). It is a highly processive enzyme with fast chain extension capability.

Platinum® Pfx DNA Polymerase is provided in inactive form, due to specific binding of the Platinum® antibody. Polymerase activity is restored after a PCR denaturation step at 94°C, providing an automatic "hot start" and increasing specificity, sensitivity, and yield (4). The high accuracy, specificity, and yield of Platinum® Pfx DNA Polymerase make it ideal for demanding PCR applications such as site-directed mutagenesis and PCR expression cloning.

For problematic and/or GC-rich templates, PCR_x Enhancer Solution is included with each kit (see the guidelines for use on page 2). The number of reactions per kit is based on a standard reaction size of 50 μl.

<u>Components</u>	<u>Kit Size</u>		
	<u>100 Rxns</u>	<u>250 Rxns</u>	<u>500 Rxns</u>
Platinum® Pfx DNA Polymerase	100 units	250 units	500 units
50 mM Magnesium Sulfate	1 ml	1 ml	1 ml
10X Pfx Amplification Buffer	1 ml	2 × 1 ml	3 × 1 ml
10X PCR _x Enhancer Solution	1 ml	2 × 1 ml	3 × 1 ml

Part no. 11708.pps

Rev. date 07/11/03

This product is distributed for laboratory research only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

For technical questions about this product, call the Invitrogen Tech-Line® U.S.A. 800 955 6288

Platinum® Pfx DNA Polymerase Storage Buffer

50 mM Tris-HCl (pH 8.0), 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, stabilizers, and 50% (v/v) glycerol.

Unit Definition

One unit of Platinum® Pfx DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 min at 74°C.

Quality Control

The DNA polymerase:antibody complex is evaluated in a DNA polymerization activity assay that measures the percent of DNA polymerase inhibition versus an uninhibited control. Platinum® Pfx DNA Polymerase is functionally tested in an amplification using 100 ng of K562 genomic DNA as a template.

Guidelines for Using PCR_x Enhancer Solution

For problematic and/or GC-rich templates, PCR_x Enhancer Solution provides higher primer specificity, a broader range of optimal magnesium concentrations, broad annealing temperatures, and improved thermostability.

Use of PCR_x Enhancer Solution is optional; use in combination with 10X Pfx Amplification Buffer, not as a substitute. PCR_x Enhancer Solution lowers the DNA melting temperature (T_m), reducing the maximum primer annealing temperature approximately 2°C per 1X PCR_x Enhancer Solution concentration, while at the same time expanding the effective annealing temperature over a much broader range. To determine the optimal reaction concentrations and conditions, we recommend starting with an annealing temperature of 55°C to 60°C and varying the amount of 10X PCR_x Enhancer Solution. For targets with higher GC content (60 to 90%), we recommend testing 10X PCR_x Enhancer Solution at final concentrations of 0.5X, 1X, 2X, and 3X.

PCR Protocol

The following procedure is suggested as a guideline and starting point when using Platinum® *Pfx* DNA Polymerase in any PCR amplification.

1. Add the following components to an autoclaved microcentrifuge tube either at ambient temperature or on ice:

<u>Component</u>	<u>Volume</u>	<u>Final Concentration</u>
10X <i>Pfx</i> Amplification Buffer	5 µl	1X
10 mM dNTP mixture*	1.5 µl	0.3 mM each
50 mM MgSO ₄	1 µl	1 mM
Primer mix (10 µM each)*	1.5 µl	0.3 µM each
Template DNA (10 pg - 200 ng)	≥1 µl	As required
Platinum® <i>Pfx</i> DNA Polymerase**	0.4–1 µl	1.0–2.5 units
Autoclaved, distilled water	to 50 µl	

*Platinum® *Pfx* DNA Polymerase will not function in reactions that contain dUTP in either the dNTP mix or the primers.

**For most targets 1 unit is sufficient. When amplifying targets above 3 kb, more enzyme may be required.

2. Mix tube contents and overlay with mineral or silicone oil, if necessary.
3. Cap the tube and centrifuge briefly to collect the contents.
4. Denature the template for 2 min at 94°C. Perform 25–35 cycles of PCR amplification as follows:

Three-step cycling

Denature: 94°C for 15 s

Anneal: 55°C for 30 s

Extend: 68°C for 1 min per kb

Two-step cycling

Denature: 94°C for 15 s

Extend: 68°C for 1 min per kb

Note: Two-step cycling can be used for long primers with high T_m.

5. Maintain the reaction at 4°C after cycling. Samples can be stored at -20°C until use.
6. Analyze the products by agarose gel electrophoresis.

References

1. Takagi, M., Nishioka, M., Kakihara, H., Kitabayashi, M., Inoue, H., Kawakami, B., Oka, M., and Imanaka, T. (1997) *Appl. Environ. Microbiol.*, 63, 4504-4510.
2. Nishioka M, Mizuguchi H, Fujiwara S, Komatsubara S, Kitabayashi M, Uemura H, Takagi M, Imanaka T. (2001) *J. Biotechnol.*, 88, 141-9.
3. Cline, J., Braman., and Hogrefe, H. H. (1996) *Nucleic Acid Res.*, 24, 3546.
4. Sharkey, D.J., Scalice, E.R., Christy, K.G., Atwood, S.M., Daiss, J.L. (1994) *BioTechnology*, 12, 506.

Limited Use Label License No. 1: Thermostable Polymerases

A license under U.S. Patents 4,683,202, 4,683,195, 4,965,188, and 5,075,216 or their foreign counterparts, owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche, Ltd. ("Roche"), has an up-front fee component and a running-royalty component. The purchase price of this product includes limited, nontransferable rights under the running-royalty component to use only this amount of the product to practice the polymerase chain reaction (PCR) and related processes described in said patents solely for the research and development activities of the purchaser when this product is used in conjunction with a thermal cycler whose use is covered by the up-front fee component. Rights to the up-front fee component must be obtained by the end user in order to have a complete license. These rights under the up-front fee component may be purchased from Applied Biosystems or obtained by purchasing an Authorized Thermal Cycler. No right to perform or offer commercial services of any kind using PCR including, without limitation, reporting the results of the purchaser's activities for a fee or other commercial consideration, is hereby granted by implication or estoppel. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

Limited Use Label License No. 14: Platinum[®] Products

Licensed to Invitrogen Corporation, under U.S. Patent Nos. 5,338,671; 5,587,287, and foreign equivalents for use in research only.

Limited Use Label License No. 37: PCR_x Products

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned by Invitrogen Corporation and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

Limited Use Label License No. 136: Pfx DNA Polymerase

This product is manufactured for Invitrogen by Toyobo Co., Ltd.

©2003 Invitrogen Corporation. All rights reserved.